

Suture Description:	Monofilament	For Smooth Passage through Tissue.
	Synthetic Absorbable	For Minimal Tissue Reaction By hydrolysis. Essentially complete in 90-120 days
Tensile strength:	Post Implantation	Post Implantation Approximate % Original
	7 days	60% (Violet) 50% (Undyed)
	14 days	30% (Violet) 20% (Undyed)
	21 days	0% (Undyed)
	28 days	0% (Violet)

Color: Violet or Undyed.





Range: 6-0 to 1 (USP). Supplied as needled sutures and ligatures.



Indications: MONOPRONE* sutures are intended for use in general soft tissue approximation and/or ligation where an absorbable material is indicated. Typical areas of use include:




- Subcuticular Sutures
- Small Intestine Anastomoses
- Urological Anastomoses



Full details in the Instruction for Use included in every package.






1/2 Circle		Taper Point Needles					Micro Needles	
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R10-180U-250	 Length: 10mm Radius: 3.18mm Curvature: 180° Chord: 6.36mm Diameter: 250µ		6-0	0.7	VIOLET	70		M6R10U70
10mm								
R12-180U-400	 Length: 12mm Radius: 3.82mm Curvature: 180° Chord: 7.64mm Diameter: 400µ		5-0	1	UNDYED	45		M5R12U45U
12mm			4-0	1.5	UNDYED	45		M4R12U45U
			4-0	1.5	UNDYED	75		M4R12U75U
R13-180U-350	 Length: 13mm Radius: 4.14mm Curvature: 180° Chord: 8.28mm Diameter: 350µ		7-0	0.5	VIOLET	45		M7R13U45
13mm			6-0	0.7	VIOLET	45	SLIM	M6R13U45S
			6-0	0.7	UNDYED	60		M6R13U60U
			6-0	0.7	VIOLET	70		M6R13U70
			6-0	0.7	UNDYED	70		M6R13U70U
			6-0	0.7	VIOLET	70	X-BLACK*	M6R13U70B
			6-0	0.7	VIOLET	70		M6R13U70
			5-0	1	VIOLET	70		M5R13U70
			5-0	1	UNDYED	70		M5R13U70U
			5-0	1	VIOLET	70		M5R13U70
			5-0	1	VIOLET	75		M5R13U75


1/2 Circle		Taper Point Needles				Micro Needles		
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R13-180U-350 13mm		Length: 13mm	Radius: 4.14mm	4-0	1.5	VIOLET	70	M4R13U70
		Curvature: 180°	Chord: 8.28mm	4-0	1.5	UNDYED	70	M4R13U70U
R13-180U-350-D 13mm		Length: 13mm	Radius: 4.14mm	6-0	0.7	VIOLET	70	M6R13U70D
		Curvature: 180°	Chord: 8.28mm	5-0	1	VIOLET	45	M5R13U45D
			Diameter: 350µ	5-0	1	VIOLET	70	M5R13U70D



3/8 Circle		Taper Point Needles				Micro Needles		
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R8-135C-230-D 8mm		Length: 8mm	Radius: 3.4mm	5-0	1	VIOLET	45	M5R8C45D
R10-135C-200-D 10mm		Length: 10mm	Radius: 4.25mm	6-0	0.7	VIOLET	90	M6R10C90D
		Curvature: 135°	Chord: 7.85mm					
R12-135C-230 12mm		Length: 12mm	Radius: 5.1mm	6-0	0.7	VIOLET	70	M6R12C70
		Curvature: 135°	Chord: 9.42mm	5-0	1	VIOLET	70	M5R12C70
			Diameter: 230µ	5-0	1	VIOLET	90	M5R12C90



3/8 Circle		Taper Point Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R12-135C-230-D 12mm		Length: 12mm	Radius: 5.1mm	VIOLET	70		M6R12C70D	
		Curvature: 135°	Chord: 9.42mm				6-0	0.7
R13-135C-530-D 13mm		Length: 13mm	Radius: 5.52mm	VIOLET	75		M5R13C75D	
		Curvature: 135°	Chord: 10.2mm				5-0	1
		Diameter: 230µ						
		Diameter: 530µ						




5/8 Circle		Taper Point Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R8-2250-300 8mm		Length: 8mm	Radius: 2.04mm	UNDYED	70		M5R8070U	
R13-2250-400 13mm		Curvature: 225°	Chord: 3.77mm				5-0	1
		Length: 13mm	Radius: 3.31mm	VIOLET	90		M4R13090	
		Diameter: 300µ						
		Diameter: 400µ						

Progressive		Taper Point Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R13-95/90P-350-B 13mm		Length: 13mm	Radius: 0mm	VIOLET	45	X-BLACK*	M5R13P45B	
		Curvature: 90/95°	Chord: 0mm				5-0	1
		Diameter: 350µ						



Progressive		Taper Cutting Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T13-95/90P-350 13mm 	Length: 13mm	5-0	1	UNDYED	45		M5T13P45U
	Radius: 0mm						
	Curvature: 95/90°						
	Chord: 0mm						
	Diameter: 350µ						



3/8 Circle		Reverse Cutting Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C12-135C-300 12mm 	Length: 12mm	6-0	0.7	VIOLET	45		M6C12C45
	Radius: 5.1mm	6-0	0.7	UNDYED	45		M6C12C45U
	Curvature: 135°	6-0	0.7	VIOLET	70		M6C12C70
	Chord: 9.42mm	5-0	1	VIOLET	45		M5C12C45
	Diameter: 300µ	5-0	1	UNDYED	60		M5C12C60U
		5-0	1	VIOLET	70		M5C12C70
		4-0	1.5	VIOLET	45		M4C12C45
C13-135C-350 13mm 	Length: 13mm	5-0	1	UNDYED	75		M5C13C75U
	Radius: 5.52mm						
	Curvature: 135°						
	Chord: 10.2mm						
	Diameter: 350µ						


Straight		Reverse Cutting Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C13-0Y-380 13mm	 Length: 13mm Curvature: 0°	4-0	1.5	UNDYED	45		M4C13Y45U
C13-0Y-380-D 13mm	 Length: 13mm Curvature: 0°	5-0	1	UNDYED	45		M5C13Y45DU
	Radius: 0mm Chord: 13mm Diameter: 380µ						



1/2 Circle		Reverse Cutting Precision Point Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PC8-180U-200 8mm	 Length: 8mm Curvature: 180°	5-0	1	VIOLET	45		M5PC8U45
		5-0	1	UNDYED	45		M5PC8U45U
		4-0	1.5	UNDYED	45		M4PC8U45U
PC9-180U-250 9mm	 Length: 9mm Curvature: 180°	6-0	0.7	UNDYED	45		M6PC9U45U
		4-0	1.5	UNDYED	45		M4PC9U45U
PC10-180U-250 10mm	 Length: 10mm Curvature: 180°	6-0	0.7	UNDYED	45		M6PC10U45U
	Radius: 2.55mm Chord: 5.1mm Diameter: 200µ						
	Radius: 2.87mm Chord: 5.74mm Diameter: 250µ						
	Radius: 3.18mm Chord: 6.36mm Diameter: 250µ						


3/8 Circle		Reverse Cutting Precision Point Needles					Micro Needles	
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
7mm		Length: 7mm Curvature: 135° Radius: 2.97mm Chord: 5.49mm Diameter: 200µ	7-0	0.5	VIOLET	45		M7PC7C45
			6-0	0.7	VIOLET	45		M6PC7C45
9mm		Length: 9mm Curvature: 135° Radius: 3.82mm Chord: 7.06mm Diameter: 240µ	3-0	2	UNDYED	45		M3PC9C45U
11mm		Length: 11mm Curvature: 135° Radius: 4.67mm Chord: 8.63mm Diameter: 350µ	7-0	0.5	VIOLET	45		M7PC11C45
			6-0	0.7	VIOLET	45		M6PC11C45
			6-0	0.7	UNDYED	45		M6PC11C45U
			5-0	1	UNDYED	25		M5PC11C25U
			5-0	1	VIOLET	45		M5PC11C45
			5-0	1	UNDYED	45		M5PC11C45U
12mm		Length: 12mm Curvature: 135° Radius: 5.1mm Chord: 9.42mm Diameter: 300µ	6-0	0.7	UNDYED	45		M6PC12C45U
			5-0	1	UNDYED	60		M5PC12C60U
13mm		Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 400µ	7-0	0.5	VIOLET	45		M7PC13C45
			6-0	0.7	VIOLET	45		M6PC13C45
			6-0	0.7	UNDYED	45		M6PC13C45U
			5-0	1	UNDYED	25		M5PC13C25U



3/8 Circle		Reverse Cutting Precision Point Needles				Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC13-135C-400 13mm		Length: 13mm	Radius: 5.52mm	5-0	1	VIOLET	45	M5PC13C45
		Curvature: 135°	Chord: 10.2mm	5-0	1	UNDYED	45	M5PC13C45U
		Diameter: 400µ	5-0	1	UNDYED	75		M5PC13C75U
		4-0	1.5	UNDYED	25		M4PC13C25U	
		4-0	1.5	VIOLET	45		M4PC13C45	
		4-0	1.5	UNDYED	45		M4PC13C45U	
		3-0	2	UNDYED	45		M3PC13C45U	
PC13-135C-400-D 13mm		Length: 13mm	Radius: 5.52mm	5-0	1	UNDYED	45	M5PC13C45DU
		Curvature: 135°	Chord: 10.2mm					
		Diameter: 400µ						



3/8 Circle		X-Prime Reverse Needles				Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X11-135C-250 11mm		Length: 11mm	Radius: 4.67mm	6-0	0.7	UNDYED	45	M6X11C45U
		Curvature: 135°	Chord: 8.63mm	5-0	1	UNDYED	45	M5X11C45U
X13-135C-350 13mm		Length: 13mm	Radius: 5.52mm	6-0	0.7	UNDYED	45	M6X13C45U
		Curvature: 135°	Chord: 10.2mm	5-0	1	VIOLET	45	M5X13C45
		Diameter: 350µ						



3/8 Circle			X-Prime Reverse Needles				Micro Needles	
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X13-135C-350 13mm		Length: 13mm Radius: 5.52mm Chord: 10.2mm Diameter: 350µ	5-0	1	UNDYED	45		M5X13C45U
			4-0	1.5	VIOLET	45		M4X13C45
			4-0	1.5	UNDYED	45		M4X13C45U



3/8 Circle			Conventional Cutting Needles				Micro Needles	
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K11-135C-300 11mm		Length: 11mm Radius: 4.67mm Chord: 8.63mm Diameter: 300µ	5-0	1	UNDYED	45		M5K11C45U
			4-0	1.5	UNDYED	45		M4K11C45U
K12-135C-300 12mm		Length: 12mm Radius: 5.1mm Chord: 9.42mm Diameter: 300µ	6-0	0.7	UNDYED	45		M6K12C45U



3/8 Circle			Conventional Cutting Precision Point Needles				Micro Needles	
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PK12-135C-300 12mm		Length: 12mm Radius: 5.1mm Chord: 9.42mm Diameter: 300µ	6-0	0.7	UNDYED	45		M6PK12C45U



3/8 Circle			X-Prime Cutting Needles					Micro Needles
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D13-135C-330 13mm		Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 330µ	6-0	0.7	UNDYED	45		M6D13C45U
			5-0	1	UNDYED	45		M5D13C45U
			4-0	1.5	UNDYED	45		M4D13C45U


1/4 Circle			Conventional Spatula Needles					Micro Needles
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Z8-112Q-330 8mm		Length: 8mm Curvature: 112° Radius: 4.09mm Chord: 6.78mm Diameter: 330µ	5-0	1	UNDYED	45		M5Z8Q45U



1/2 Circle			X-Trocar Needles					Micro Needles
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E13-180U-300 13mm		Length: 13mm Curvature: 180° Radius: 4.14mm Chord: 8.28mm Diameter: 300µ	6-0	0.7	VIOLET	45		M6E13U45
			6-0	0.7	VIOLET	70		M6E13U70
			5-0	1	VIOLET	70		M5E13U70




1/2 Circle		Taper Point Plus Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Q17U		6-0	0.7	VIOLET	70		M6Q17U70
		5-0	1	VIOLET	70		M5Q17U70
		4-0	1.5	VIOLET	70		M4Q17U70
		4-0	1.5	VIOLET	70	X-BLACK*	M4Q17U70B
		3-0	2	VIOLET	70		M3Q17U70
		3-0	2	UNDYED	70		M3Q17U70U
Q17U-D-B		5-0	1	VIOLET	70	X-BLACK*	M5Q17U70DB
		5-0	1	VIOLET	70	X-BLACK*	M5Q17U70DB
Q22U		4-0	1.5	VIOLET	70		M4Q22U70
		4-0	1.5	VIOLET	70	X-BLACK*	M4Q22U70B
		3-0	2	VIOLET	70		M3Q22U70
		3-0	2	VIOLET	70	X-BLACK*	M3Q22U70B
		2-0	3	VIOLET	70	X-BLACK*	M2Q22U70B
Q22U-D-B		4-0	1.5	VIOLET	70	X-BLACK*	M4Q22U70DB
		4-0	1.5	VIOLET	90	X-BLACK*	M4Q22U90DB
		3-0	2	VIOLET	70	X-BLACK*	M3Q22U70DB



1/2 Circle		Taper Point Plus Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q26U  26mm	4-0	1.5	VIOLET	70		M4Q26U70	
	4-0	1.5	VIOLET	70	X-BLACK*	M4Q26U70B	
	3-0	2	VIOLET	70		M3Q26U70	
	3-0	2	UNDYED	70		M3Q26U70U	
	3-0	2	VIOLET	70	X-BLACK*	M3Q26U70B	
	3-0	2	VIOLET	8 X 45	"CONTROLRELEASE",X-BLACK**	M3Q26U70B	
	3-0	2	VIOLET	8 X 70	"CONTROLRELEASE",X-BLACK**	M3Q26U70B	
	2-0	3	VIOLET	70		M2Q26U70	
	2-0	3	VIOLET	70	X-BLACK*	M2Q26U70B	
	2-0	3	VIOLET	70	HEAVY NEEDLE	M2Q26U70H	
	0	3.5	VIOLET	70		M0Q26U70	
	0	3.5	VIOLET	70	HEAVY NEEDLE	M0Q26U70H	
Q26U-D-B  26mm	3-0	2	VIOLET	70	X-BLACK*	M3Q26U70DB	
	3-0	2	VIOLET	90	X-BLACK*	M3Q26U90DB	
	2-0	3	VIOLET	70	X-BLACK*	M2Q26U70DB	


1/2 Circle		Taper Point Plus Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Q31U		31mm	3-0	2	VIOLET	70	M3Q31U70
			2-0	3	VIOLET	70	M2Q31U70
Q36U		36mm	2-0	3	VIOLET	70	M2Q36U70
			2-0	3	VIOLET	70	HEAVY NEEDLE M2Q36U70H
			0	3.5	VIOLET	70	HEAVY NEEDLE M0Q36U70H
			0	3.5	VIOLET	90	HEAVY NEEDLE M0Q36U90H
			1	4	VIOLET	70	HEAVY NEEDLE M1Q36U70H



1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R17U		17mm	6-0	0.7	UNDYED	70	M6R17U70U
			6-0	0.7	VIOLET	75	M6R17U75
			5-0	1	VIOLET	70	X-BLACK* M5R17U70B
			5-0	1	VIOLET	75	M5R17U75
			5-0	1	UNDYED	75	M5R17U75U
			4-0	1.5	VIOLET	45	M4R17U45



1/2 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R17U  17mm	4-0	1.5	VIOLET	75		M4R17U75	
	4-0	1.5	UNDYED	75		M4R17U75U	
	4-0	1.5	VIOLET	90		M4R17U90	
	3-0	2	VIOLET	45		M3R17U45	
	3-0	2	VIOLET	70	X-BLACK*	M3R17U70B	
	3-0	2	VIOLET	75		M3R17U75	
	3-0	2	UNDYED	75		M3R17U75U	
	2-0	3	VIOLET	45	X-BLACK*	M2R17U45B	
	2-0	3	VIOLET	75		M2R17U75	
R17U-D  17mm	5-0	1	VIOLET	70		M5R17U70D	
	4-0	1.5	VIOLET	70		M4R17U70D	
	4-0	1.5	VIOLET	70	X-BLACK*	M4R17U70DB	
	4-0	1.5	VIOLET	75		M4R17U75D	
	4-0	1.5	VIOLET	90		M4R17U90D	
	3-0	2	VIOLET	70		M3R17U70D	
	3-0	2	VIOLET	75		M3R17U75D	



1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R20U		20mm	4-0	1.5	UNDYED	75	M4R20U75U
R22U		22mm	6-0	0.7	VIOLET	45	M6R22U45
			5-0	1	VIOLET	45	M5R22U45
			5-0	1	VIOLET	70	M5R22U70
			4-0	1.5	VIOLET	75	M4R22U75
			4-0	1.5	UNDYED	75	M4R22U75U
			3-0	2	VIOLET	30	M3R22U30
			3-0	2	UNDYED	45	M3R22U45U
			3-0	2	VIOLET	70	M3R22U70
			3-0	2	UNDYED	70	M3R22U70U
			3-0	2	VIOLET	120	M3R22U120
			2-0	3	VIOLET	30	M2R22U30
			2-0	3	UNDYED	75	M2R22U75U
			2-0	3	VIOLET	120	M2R22U120
			R22U-D		22mm	4-0	1.5
4-0	1.5	VIOLET				75	M4R22U75D



1/2 Circle		Taper Point Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R22U-D		22mm	4-0	1.5	VIOLET	90		M4R22U90D
			3-0	2	VIOLET	70		M3R22U70D
			3-0	2	VIOLET	90		M3R22U90D
R26U		26mm	5-0	1	VIOLET	90		M5R26U90
			4-0	1.5	VIOLET	45		M4R26U45
			4-0	1.5	VIOLET	70	HEAVY NEEDLE	M4R26U70H
			4-0	1.5	VIOLET	75		M4R26U75
			4-0	1.5	UNDYED	75		M4R26U75U
			4-0	1.5	VIOLET	90		M4R26U90
			3-0	2	VIOLET	30		M3R26U30
			3-0	2	VIOLET	45		M3R26U45
			3-0	2	VIOLET	45	X-BLACK*	M3R26U45B
			3-0	2	VIOLET	70		M3R26U70
			3-0	2	VIOLET	70	HEAVY NEEDLE	M3R26U70H
			3-0	2	UNDYED	75		M3R26U75U
			3-0	2	VIOLET	90		M3R26U90



1/2 Circle		Taper Point Needles				
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R26U  26mm	3-0	2	VIOLET	90	HEAVY NEEDLE	M3R26U90H
	3-0	2	VIOLET	120		M3R26U120
	2-0	3	VIOLET	30		M2R26U30
	2-0	3	VIOLET	45		M2R26U45
	2-0	3	VIOLET	75		M2R26U75
	2-0	3	UNDYED	75		M2R26U75U
	2-0	3	VIOLET	90		M2R26U90
	2-0	3	VIOLET	90	HEAVY NEEDLE	M2R26U90H
	2-0	3	VIOLET	120		M2R26U120
	0	3.5	UNDYED	45		M0R26U45U
	0	3.5	VIOLET	70		M0R26U70
	0	3.5	UNDYED	75		M0R26U75U
	0	3.5	VIOLET	90	HEAVY NEEDLE	M0R26U90H
	1	4	VIOLET	70		M1R26U70
	1	4	VIOLET	70	HEAVY NEEDLE	M1R26U70H
	1	4	VIOLET	90	HEAVY NEEDLE	M1R26U90H


1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R26U		26mm	2	5	VIOLET	70	HEAVY NEEDLE	MXR26U70H
			R26U-D		26mm	4-0	1.5	VIOLET
4-0	1.5	VIOLET	90			M4R26U90D		
3-0	2	VIOLET	70			M3R26U70D		
3-0	2	VIOLET	90			M3R26U90D		
2-0	3	VIOLET	20			M2R26U20D		
2-0	3	VIOLET	90			M2R26U90D		
R27U		27mm	3-0	2	VIOLET	75	M3R27U75	
			2-0	3	VIOLET	75	M2R27U75	
			2-0	3	UNDYED	75	M2R27U75U	
			0	3.5	VIOLET	75	M0R27U75	
			0	3.5	UNDYED	75	M0R27U75U	
			1	4	VIOLET	75	M1R27U75	
			1	4	UNDYED	75	M1R27U75U	




1/2 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R30U  30mm	4-0	1.5	VIOLET	70		M4R30U70	
	3-0	2	VIOLET	45		M3R30U45	
	3-0	2	VIOLET	75		M3R30U75	
	3-0	2	UNDYED	75		M3R30U75U	
	2-0	3	VIOLET	75		M2R30U75	
	2-0	3	VIOLET	90		M2R30U90	
	0	3.5	VIOLET	70		M0R30U70	
	0	3.5	VIOLET	70	HEAVY NEEDLE	M0R30U70H	
	0	3.5	VIOLET	90		M0R30U90	
	1	4	VIOLET	70		M1R30U70	
	1	4	VIOLET	70	HEAVY NEEDLE	M1R30U70H	
R30U-D  30mm	3-0	2	VIOLET	70		M3R30U70D	
	3-0	2	VIOLET	90		M3R30U90D	
	2-0	3	VIOLET	70		M2R30U70D	
	2-0	3	VIOLET	90		M2R30U90D	



1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R31U		2-0	3	UNDYED	75		M2R31U75U
		0	3.5	VIOLET	45		M0R31U45
		0	3.5	UNDYED	75		M0R31U75U
		1	4	UNDYED	75		M1R31U75U
R36U		3-0	2	VIOLET	45		M3R36U45
		3-0	2	VIOLET	70		M3R36U70
		3-0	2	VIOLET	90		M3R36U90
		3-0	2	UNDYED	90		M3R36U90U
		2-0	3	UNDYED	75		M2R36U75U
		0	3.5	VIOLET	45		M0R36U45
		0	3.5	UNDYED	45		M0R36U45U
		0	3.5	VIOLET	70		M0R36U70
		0	3.5	VIOLET	90		M0R36U90
		0	3.5	UNDYED	90		M0R36U90U
		1	4	VIOLET	70		M1R36U70
		1	4	VIOLET	90		M1R36U90




1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R36U		36mm	1	4	UNDYED	90		M1R36U90U
			R37U		37mm	4-0	1.5	VIOLET
3-0	2	VIOLET	70				M3R37U70	
3-0	2	VIOLET	70			HEAVY NEEDLE	M3R37U70H	
3-0	2	VIOLET	90				M3R37U90	
3-0	2	UNDYED	90				M3R37U90U	
3-0	2	VIOLET	90			HEAVY NEEDLE	M3R37U90H	
2-0	3	VIOLET	45				M2R37U45	
2-0	3	VIOLET	70			HEAVY NEEDLE	M2R37U70H	
2-0	3	VIOLET	75				M2R37U75	
2-0	3	UNDYED	75				M2R37U75U	
2-0	3	VIOLET	90				M2R37U90	
2-0	3	UNDYED	90				M2R37U90U	
2-0	3	VIOLET	90			HEAVY NEEDLE	M2R37U90H	
0	3.5	VIOLET	45				M0R37U45	



1/2 Circle		Taper Point Needles			◎			
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R37U		37mm	0	3.5	UNDYED	75		M0R37U75U
			0	3.5	VIOLET	75	HEAVY NEEDLE	M0R37U75H
			0	3.5	VIOLET	90		M0R37U90
			0	3.5	UNDYED	90		M0R37U90U
			0	3.5	VIOLET	90	HEAVY NEEDLE	M0R37U90H
			1	4	VIOLET	70	HEAVY NEEDLE	M1R37U70H
			1	4	VIOLET	75		M1R37U75
			1	4	UNDYED	75		M1R37U75U
			1	4	VIOLET	90		M1R37U90
			1	4	UNDYED	90		M1R37U90U
			1	4	VIOLET	90	HEAVY NEEDLE	M1R37U90H
			1	4	UNDYED	120	HEAVY NEEDLE	M1R37U120HU
			2	5	VIOLET	70		MXR37U70
R40U		40mm	3-0	2	VIOLET	70		M3R40U70
			3-0	2	UNDYED	90		M3R40U90U
			2-0	3	VIOLET	70	HEAVY NEEDLE	M2R40U70H





1/2 Circle		Taper Point Needles				
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R40U  40mm	2-0	3	VIOLET	75		M2R40U75
	2-0	3	UNDYED	75		M2R40U75U
	2-0	3	UNDYED	90		M2R40U90U
	2-0	3	VIOLET	90	HEAVY NEEDLE	M2R40U90H
	0	3.5	VIOLET	45		M0R40U45
	0	3.5	VIOLET	45	HEAVY NEEDLE	M0R40U45H
	0	3.5	VIOLET	70	HEAVY NEEDLE	M0R40U70H
	0	3.5	UNDYED	75		M0R40U75U
	0	3.5	VIOLET	90		M0R40U90
	0	3.5	UNDYED	90		M0R40U90U
	0	3.5	VIOLET	90	HEAVY NEEDLE	M0R40U90H
	0	3.5	UNDYED	90	HEAVY NEEDLE	M0R40U90HU
	0	3.5	VIOLET	150	LOOP	M0R40U150L
	1	4	VIOLET	70		M1R40U70
	1	4	VIOLET	70	HEAVY NEEDLE	M1R40U70H
	1	4	UNDYED	75		M1R40U75U



1/2 Circle		Taper Point Needles			◎		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R40U		1	4	VIOLET	90		M1R40U90
		1	4	UNDYED	90		M1R40U90U
		1	4	VIOLET	90	HEAVY NEEDLE	M1R40U90H
		1	4	VIOLET	150	LOOP	M1R40U150L
		2	5	VIOLET	70	HEAVY NEEDLE	MXR40U70H
		2	5	VIOLET	90	HEAVY NEEDLE	MXR40U90H
R43U		0	3.5	VIOLET	70		M0R43U70
		0	3.5	VIOLET	70	HEAVY NEEDLE	M0R43U70H
		1	4	VIOLET	90	HEAVY NEEDLE	M1R43U90H
R48U		3-0	2	UNDYED	90		M3R48U90U
		2-0	3	UNDYED	75		M2R48U75U
		2-0	3	VIOLET	90		M2R48U90
		2-0	3	VIOLET	90	HEAVY NEEDLE	M2R48U90H
		0	3.5	VIOLET	90		M0R48U90
		0	3.5	UNDYED	90		M0R48U90U
		0	3.5	UNDYED	120		M0R48U120U




1/2 Circle		Taper Point Needles			◎		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R48U  48mm		0	3.5	VIOLET	150	LOOP	M0R48U150L
		0	3.5	VIOLET	150		M0R48U150
		1	4	VIOLET	45		M1R48U45
		1	4	UNDYED	70		M1R48U70U
		1	4	VIOLET	90		M1R48U90
		1	4	UNDYED	90		M1R48U90U
		1	4	VIOLET	90	HEAVY NEEDLE	M1R48U90H
		1	4	VIOLET	120		M1R48U120
		1	4	VIOLET	140		M1R48U140
		1	4	VIOLET	150	LOOP	M1R48U150L
		2	5	VIOLET	90		MXR48U90
R60U  60mm		1	4	VIOLET	90		M1R60U90





1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R65U		65mm	2-0	3	VIOLET	90		M2R65U90
			0	3.5	VIOLET	150	LOOP	M0R65U150L
			1	4	VIOLET	150	LOOP	M1R65U150L
R70U		70mm	3-0	2	UNDYED	90		M3R70U90U
			1	4	VIOLET	90		M1R70U90
R76U		76mm	1	4	VIOLET	90		M1R76U90
			2	5	VIOLET	150		MXR76U150





3/8 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R17C		17mm	4-0	1.5	UNDYED	70		M4R17C70U
R18C		18mm	5-0	1	VIOLET	70		M5R18C70
			4-0	1.5	VIOLET	70		M4R18C70
			3-0	2	VIOLET	30	LOOP	M3R18C30L



3/8 Circle		Taper Point Needles			◎		
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R18C-D	 18mm	6-0	0.7	VIOLET	70		M6R18C70D
R22C	 22mm	4-0	1.5	VIOLET	70		M4R22C70
		4-0	1.5	VIOLET	120		M4R22C120
		3-0	2	VIOLET	70		M3R22C70
		3-0	2	VIOLET	120		M3R22C120
		2-0	3	VIOLET	30		M2R22C30
		2-0	3	VIOLET	120		M2R22C120
R26C	 26mm	3-0	2	VIOLET	70		M3R26C70
		3-0	2	VIOLET	120		M3R26C120
		2-0	3	VIOLET	70		M2R26C70
R26C-D	 26mm	3-0	2	VIOLET	70		M3R26C70D
R32C	 32mm	4-0	1.5	VIOLET	70		M4R32C70
		3-0	2	VIOLET	70		M3R32C70
		2-0	3	VIOLET	70		M2R32C70


3/8 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R32C		32mm	0	3.5	VIOLET	70	M0R32C70
			1	4	VIOLET	70	M1R32C70
R32C-D		32mm	2-0	3	VIOLET	70	M2R32C70D





5/8 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R170		17mm	5-0	1	VIOLET	90	M5R17090
			4-0	1.5	VIOLET	90	M4R17090
R220		22mm	5-0	1	VIOLET	90	M5R22090
			4-0	1.5	VIOLET	90	M4R22090
R260		26mm	4-0	1.5	VIOLET	70	M4R26070
			3-0	2	VIOLET	70	M3R26070
			3-0	2	UNDYED	70	M3R26070U
			3-0	2	VIOLET	120	M3R260120
			2-0	3	VIOLET	70	M2R26070


5/8 Circle		Taper Point Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R260		26mm	2-0	3	UNDYED	70	M2R26070U
			2-0	3	VIOLET	90	M2R26090
			0	3.5	VIOLET	70	M0R26070
R260-D		26mm	4-0	1.5	VIOLET	70	M4R26070D
			3-0	2	VIOLET	70	M3R26070D
			3-0	2	VIOLET	90	M3R26090D
			2-0	3	VIOLET	45	M2R26045D
			2-0	3	VIOLET	70	M2R26070D
			2-0	3	VIOLET	90	M2R26090D
R270		27mm	3-0	2	VIOLET	75	M3R27075
			2-0	3	UNDYED	75	M2R27075U
			0	3.5	VIOLET	75	M0R27075
			0	3.5	UNDYED	75	M0R27075U
R360		36mm	3-0	2	VIOLET	70	M3R36070
			2-0	3	VIOLET	70	M2R36070
			0	3.5	VIOLET	70	M0R36070







5/8 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R360-D		36mm	3-0	2	VIOLET	90	M3R36090D
			2-0	3	UNDYED	70	M2R36070DU
			2-0	3	VIOLET	90	M2R36090D
			0	3.5	VIOLET	90	M0R36090D
R370		37mm	2-0	3	VIOLET	75	M2R37075
			2-0	3	UNDYED	75	M2R37075U
R370-D		37mm	3-0	2	VIOLET	75	M3R37075D
R400			2-0	3	VIOLET	70	M2R40070
			0	3.5	VIOLET	75	M0R40075
		40mm	0	3.5	VIOLET	90	M0R40090


Ski		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R19S		19mm	3-0	2	VIOLET	120	M3R19S120
R26S		26mm	3-0	2	VIOLET	30	M3R26S30


Ski		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R26S  26mm		3-0	2	VIOLET	120		M3R26S120
		2-0	3	VIOLET	30		M2R26S30
		2-0	3	VIOLET	120		M2R26S120



Straight		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R19Y  19mm		4-0	1.5	VIOLET	30		M4R19Y30
		3-0	2	VIOLET	30		M3R19Y30
R22Y  22mm		0	3.5	UNDYED	90		M0R22Y90U
R65Y  65mm		3-0	2	VIOLET	90		M3R65Y90
R65Y-D  65mm		2-0	3	VIOLET	70		M2R65Y70D



1/2 Circle		X-GUARD Blunt Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
G70U  70mm		1	4	VIOLET	90		M1G70U90





1/2 Circle		Blunt Point Needles			○			
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
B26U	 26mm		3-0	2	UNDYED	75		M3B26U75U
B35U	 35mm		2-0	3	VIOLET	70		M2B35U70
			0	3.5	VIOLET	70		MOB35U70
			0	3.5	VIOLET	90		MOB35U90
B36U	 36mm		3-0	2	UNDYED	90		M3B36U90U
			0	3.5	VIOLET	70		MOB36U70
			0	3.5	UNDYED	90		MOB36U90U
			1	4	UNDYED	90		M1B36U90U
B37U	 37mm		2-0	3	VIOLET	75		M2B37U75
			0	3.5	VIOLET	75		MOB37U75
			0	3.5	UNDYED	90		MOB37U90U
B40U	 40mm		1	4	VIOLET	90		M1B40U90
B48U	 48mm		0	3.5	UNDYED	90		MOB48U90U
			0	3.5	VIOLET	150	LOOP	MOB48U150L
			1	4	UNDYED	90		M1B48U90U





1/2 Circle		Blunt Point Needles			○			
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
B50U		50mm	0	3.5	UNDYED	90		M0B50U90U
			1	4	VIOLET	90		M1B50U90




3/8 Circle		Blunt Point Needles			○			
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
B64C		64mm	1	4	VIOLET	90		M1B64C90



1/2 Circle		Taper Cutting Needles			◕			
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
T26U		26mm	2-0	3	VIOLET	70		M2T26U70
T36U		36mm	2-0	3	VIOLET	70		M2T36U70
			0	3.5	VIOLET	70		M0T36U70


1/2 Circle		Taper Cutting Needles			▼		
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T36U	 36mm	1	4	VIOLET	90		M1T36U90
T37U	 37mm	2-0	3	UNDYED	90		M2T37U90U
		1	4	UNDYED	90		M1T37U90U



1/2 Circle		Reverse Cutting Needles			▼		
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C18U	 18mm	5-0	1	UNDYED	70		M5C18U70U
		4-0	1.5	VIOLET	70		M4C18U70
C18U-D	 18mm	4-0	1.5	UNDYED	30		M4C18U30DU
C22U	 22mm	3-0	2	UNDYED	90		M3C22U90U
		2-0	3	UNDYED	70		M2C22U70U
C23U	 23mm	4-0	1.5	UNDYED	70		M4C23U70U
		3-0	2	UNDYED	70		M3C23U70U
		2-0	3	UNDYED	70		M2C23U70U



1/2 Circle		Reverse Cutting Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C24U		24mm	4-0	1.5	VIOLET	75		M4C24U75
			3-0	2	VIOLET	75		M3C24U75
C26U		26mm	3-0	2	VIOLET	70		M3C26U70
			3-0	2	VIOLET	70	HEAVY NEEDLE	M3C26U70H
			2-0	3	VIOLET	70		M2C26U70
			2-0	3	UNDYED	70		M2C26U70U
			2-0	3	VIOLET	70	HEAVY NEEDLE	M2C26U70H
			2-0	3	VIOLET	90		M2C26U90
			0	3.5	VIOLET	70	HEAVY NEEDLE	M0C26U70H
			0	3.5	UNDYED	70	HEAVY NEEDLE	M0C26U70HU
C30U		30mm	2-0	3	UNDYED	70		M2C30U70U
			0	3.5	VIOLET	70		M0C30U70
C36U		36mm	2-0	3	UNDYED	70		M2C36U70U
			0	3.5	VIOLET	70		M0C36U70
			0	3.5	UNDYED	70		M0C36U70U
			0	3.5	VIOLET	90		M0C36U90






1/2 Circle		Reverse Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C37U 	37mm	3-0	2	UNDYED	70		M3C37U70U
		3-0	2	VIOLET	90		M3C37U90
		2-0	3	VIOLET	70	HEAVY NEEDLE	M2C37U70H
		2-0	3	VIOLET	75		M2C37U75
		2-0	3	UNDYED	75		M2C37U75U
		2-0	3	VIOLET	90	HEAVY NEEDLE	M2C37U90H
		0	3.5	VIOLET	70	HEAVY NEEDLE	M0C37U70H
		0	3.5	VIOLET	75		M0C37U75
		0	3.5	UNDYED	75		M0C37U75U
		0	3.5	VIOLET	90	HEAVY NEEDLE	M0C37U90H
		1	4	VIOLET	70	HEAVY NEEDLE	M1C37U70H
		1	4	VIOLET	75		M1C37U75
C37U-D 	37mm	3-0	2	UNDYED	135		M3C37U135DU
		2-0	3	UNDYED	150	LOOP	M2C37U150DUL
C40U 	40mm	0	3.5	UNDYED	90		M0C40U90U
		1	4	UNDYED	90		M1C40U90U




1/2 Circle		Reverse Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C40U  40mm	2	5	VIOLET	90		MXC40U90	
	0	3.5	UNDYED	90		M0C48U90U	
C48U  48mm	1	4	VIOLET	90		M1C48U90	
	1	4	UNDYED	90		M1C48U90U	
	2	5	VIOLET	90		MXC48U90	


3/8 Circle		Reverse Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C16C  16mm	6-0	0.7	UNDYED	45		M6C16C45U	
	5-0	1	VIOLET	45		M5C16C45	
	5-0	1	UNDYED	70		M5C16C70U	
	5-0	1	UNDYED	75		M5C16C75U	
	4-0	1.5	VIOLET	45		M4C16C45	
	4-0	1.5	UNDYED	75		M4C16C75U	



3/8 Circle		Reverse Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C16C		16mm	3-0	2	VIOLET	45	M3C16C45
			3-0	2	UNDYED	70	M3C16C70U
C19C		19mm	5-0	1	VIOLET	45	M5C19C45
			5-0	1	UNDYED	45	M5C19C45U
			5-0	1	UNDYED	75	M5C19C75U
			4-0	1.5	VIOLET	45	M4C19C45
			4-0	1.5	UNDYED	45	M4C19C45U
			4-0	1.5	VIOLET	75	M4C19C75
			4-0	1.5	UNDYED	75	M4C19C75U
			4-0	1.5	VIOLET	90	M4C19C90
			3-0	2	VIOLET	45	M3C19C45
			3-0	2	UNDYED	45	M3C19C45U
			3-0	2	UNDYED	75	M3C19C75U
			3-0	2	VIOLET	90	M3C19C90
			2-0	3	UNDYED	45	M2C19C45U
			2-0	3	VIOLET	70	M2C19C70


3/8 Circle		Reverse Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C24C		24mm	4-0	1.5	UNDYED	45		M4C24C45U
			4-0	1.5	VIOLET	70		M4C24C70
			4-0	1.5	UNDYED	70		M4C24C70U
			3-0	2	VIOLET	45		M3C24C45
			3-0	2	VIOLET	70		M3C24C70
			3-0	2	UNDYED	70		M3C24C70U
			3-0	2	UNDYED	75		M3C24C75U
			3-0	2	VIOLET	90		M3C24C90
			2-0	3	UNDYED	45		M2C24C45U
			2-0	3	UNDYED	70		M2C24C70U
			2-0	3	VIOLET	75		M2C24C75
			2-0	3	VIOLET	90		M2C24C90
			0	3.5	VIOLET	70		M0C24C70
			1	4	UNDYED	70		M1C24C70U
C26C		26mm	4-0	1.5	UNDYED	75		M4C26C75U
			3-0	2	UNDYED	75		M3C26C75U


3/8 Circle		Reverse Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C26C	 26mm	2-0	3	UNDYED	75		M2C26C75U
C30C	 30mm	4-0	1.5	UNDYED	70		M4C30C70U
		3-0	2	UNDYED	75		M3C30C75U
		2-0	3	UNDYED	45		M2C30C45U
		2-0	3	VIOLET	70		M2C30C70
		0	3.5	VIOLET	70		M0C30C70
		0	3.5	UNDYED	75		M0C30C75U
		0	3.5	VIOLET	90		M0C30C90
		1	4	VIOLET	70		M1C30C70
C31C	 31mm	2-0	3	UNDYED	75		M2C31C75U
C36C	 36mm	2-0	3	UNDYED	75		M2C36C75U
C37C	 37mm	0	3.5	UNDYED	75		M0C37C75U
		0	3.5	UNDYED	75	HEAVY NEEDLE	M0C37C75HU





3/8 Circle		Reverse Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C37C	 37mm	1	4	UNDYED	75		M1C37C75U
C39C	 39mm	2-0	3	UNDYED	70		M2C39C70U
C40C	 40mm	0	3.5	UNDYED	75		M0C40C75U
		0	3.5	UNDYED	75	HEAVY NEEDLE	M0C40C75HU
		1	4	UNDYED	75		M1C40C75U
		1	4	UNDYED	75	HEAVY NEEDLE	M1C40C75HU







5/8 Circle		Reverse Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C250	 25mm	3-0	2	UNDYED	70		M3C25070U



Straight		Reverse Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C51Y  51mm		4-0	1.5	UNDYED	70		M4C51Y70U
		3-0	2	UNDYED	70		M3C51Y70U
		2-0	3	UNDYED	70		M2C51Y70U
C60Y  60mm		4-0	1.5	UNDYED	70		M4C60Y70U
		3-0	2	VIOLET	70		M3C60Y70
		2-0	3	UNDYED	70		M2C60Y70U

1/2 Circle		Reverse Cutting Precision Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PC15U  15mm		4-0	1.5	UNDYED	45		M4PC15U45U



3/8 Circle		Reverse Cutting Precision Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PC16C  16mm		5-0	1	VIOLET	45		M5PC16C45
		5-0	1	UNDYED	70		M5PC16C70U


3/8 Circle		Reverse Cutting Precision Point Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC16C		16mm	4-0	1.5	UNDYED	45		M4PC16C45U
			4-0	1.5	VIOLET	70		M4PC16C70
PC19C		19mm	5-0	1	UNDYED	25		M5PC19C25U
			5-0	1	UNDYED	45		M5PC19C45U
			5-0	1	UNDYED	75		M5PC19C75U
			4-0	1.5	UNDYED	25		M4PC19C25U
			4-0	1.5	VIOLET	45		M4PC19C45
			4-0	1.5	UNDYED	45		M4PC19C45U
			4-0	1.5	VIOLET	70		M4PC19C70
			4-0	1.5	UNDYED	70		M4PC19C70U
			3-0	2	UNDYED	25		M3PC19C25U
			3-0	2	UNDYED	45		M3PC19C45U
			3-0	2	UNDYED	75		M3PC19C70U
PC24C		24mm	4-0	1.5	UNDYED	45		M4PC24C45U
			4-0	1.5	UNDYED	70		M4PC24C70U
			3-0	2	UNDYED	45		M3PC24C45U



3/8 Circle		Reverse Cutting Precision Point Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PC24C		3-0	2	VIOLET	70		M3PC24C70
		3-0	2	UNDYED	70	HEAVY NEEDLE	M3PC24C70HU
		2-0	3	UNDYED	45		M2PC24C45U
		2-0	3	UNDYED	70		M2PC24C70U
PC26C		4-0	1.5	UNDYED	75		M4PC26C75U
		3-0	2	UNDYED	75		M3PC26C75U
		2-0	3	UNDYED	75		M2PC26C75U
PC30C		3-0	2	UNDYED	75		M3PC30C75U
PC31C		2-0	3	UNDYED	75		M2PC31C75U
PC36C		2-0	3	UNDYED	75		M2PC36C75U






1/2 Circle		X-Prime Reverse Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X16U	 16mm	4-0	1.5	UNDYED	45		M4X16U45U



3/8 Circle		X-Prime Reverse Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X16C	 16mm	5-0	1	UNDYED	45		M5X16C45U
		4-0	1.5	UNDYED	45		M4X16C45U
X19C	 19mm	5-0	1	UNDYED	45		M5X19C45U
		4-0	1.5	VIOLET	45		M4X19C45
		4-0	1.5	UNDYED	45		M4X19C45U
		4-0	1.5	UNDYED	70		M4X19C70U
		3-0	2	UNDYED	45		M3X19C45U
		3-0	2	UNDYED	70		M3X19C70U
		2-0	3	UNDYED	70		M2X19C70U
X24C	 24mm	4-0	1.5	UNDYED	45		M4X24C45U
		4-0	1.5	UNDYED	70		M4X24C70U
		4-0	1.5	UNDYED	90		M4X24C90U


3/8 Circle		X-Prime Reverse Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X24C		24mm	3-0	2	UNDYED	70	M3X24C70U
			2-0	3	UNDYED	70	M2X24C70U
X26C		26mm	4-0	1.5	UNDYED	70	M4X26C70U
			3-0	2	UNDYED	70	M3X26C70U




Progressive		X-Prime Reverse Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X16P		16mm	5-0	1	UNDYED	45	M5X16P45U
			4-0	1.5	UNDYED	90	M4X16P90U




3/8 Circle		Conventional Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K16C		16mm	5-0	1	UNDYED	75	M5K16C75U
			4-0	1.5	UNDYED	75	M4K16C75U
K19C		19mm	5-0	1	UNDYED	75	M5K19C75U
			4-0	1.5	UNDYED	75	M4K19C75U



3/8 Circle		Conventional Cutting Needles						
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K19C	 19mm		3-0	2	UNDYED	75		M3K19C75U
K26C	 26mm		4-0	1.5	UNDYED	75		M4K26C75U
			3-0	2	UNDYED	75		M3K26C75U
			2-0	3	UNDYED	75		M2K26C75U
K30C	 30mm		3-0	2	UNDYED	75		M3K30C75U
K31C	 31mm		2-0	3	UNDYED	75		M2K31C75U
K36C	 36mm		2-0	3	UNDYED	75		M2K36C75U




Straight		Conventional Cutting Needles						
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K48Y	 48mm		4-0	1.5	UNDYED	70		M4K48Y70U
K51Y	 51mm		3-0	2	UNDYED	75		M3K51Y75U
			2-0	3	UNDYED	75		M2K51Y75U



Straight		Conventional Cutting Needles			▲		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K60Y	 60mm	4-0	1.5	UNDYED	70		M4K60Y70U
		3-0	2	UNDYED	70		M3K60Y70U







3/8 Circle		Conventional Cutting Precision Point Needles			▲		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PK16C	 16mm	5-0	1	UNDYED	25		M5PK16C25U
		5-0	1	UNDYED	45		M5PK16C45U
		5-0	1	UNDYED	75		M5PK16C75U
		4-0	1.5	UNDYED	25		M4PK16C25U
		4-0	1.5	UNDYED	45		M4PK16C45U
		4-0	1.5	UNDYED	75		M4PK16C75U
PK19C	 19mm	5-0	1	UNDYED	75		M5PK19C75U
		4-0	1.5	UNDYED	75		M4PK19C75U
		3-0	2	UNDYED	75		M3PK19C75U
PK26C	 26mm	4-0	1.5	UNDYED	75		M4PK26C75U
		3-0	2	UNDYED	75		M3PK26C75U
		2-0	3	UNDYED	75		M2PK26C75U






3/8 Circle		Conventional Cutting Precision Point Needles						
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PK31C	 31mm		2-0	3	UNDYED	75		M2PK31C75U
PK36C	 36mm		2-0	3	UNDYED	75		M2PK36C75U





1/2 Circle		X-Prime Cutting Needles						
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D16U	 16mm		4-0	1.5	UNDYED	45		M4D16U45U






3/8 Circle		X-Prime Cutting Needles						
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D16C	 16mm		5-0	1	UNDYED	45		M5D16C45U
			4-0	1.5	UNDYED	45		M4D16C45U
D19C	 19mm		5-0	1	UNDYED	45		M5D19C45U
			4-0	1.5	UNDYED	45		M4D19C45U



5/8 Circle		X-Prime Cutting Needles						
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D160	 16mm		4-0	1.5	UNDYED	45		M4D16045U


1/2 Circle		X-Trocar Needles						
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E17U		17mm	6-0	0.7	VIOLET	70		M6E17U70
			5-0	1	VIOLET	70		M5E17U70
			4-0	1.5	VIOLET	70		M4E17U70
			3-0	2	VIOLET	70		M3E17U70
			2-0	3	VIOLET	70		M2E17U70
E17U-D	 	17mm	4-0	1.5	VIOLET	70		M4E17U70D
			3-0	2	VIOLET	70		M3E17U70D
E22U-D	 	22mm	4-0	1.5	VIOLET	70		M4E22U70D
			4-0	1.5	VIOLET	90		M4E22U90D
			3-0	2	VIOLET	70		M3E22U70D
			3-0	2	VIOLET	90		M3E22U90D


1/2 Circle		X-Trocar Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
E26U		26mm	4-0	1.5	VIOLET	70		M4E26U70
			3-0	2	VIOLET	70		M3E26U70
			2-0	3	VIOLET	70		M2E26U70
			0	3.5	VIOLET	70		M0E26U70
			0	3.5	VIOLET	90		M0E26U90
			1	4	VIOLET	70		M1E26U70
			1	4	VIOLET	90		M1E26U90
E26U-D		26mm	4-0	1.5	VIOLET	70		M4E26U70D
			3-0	2	VIOLET	70		M3E26U70D
E30U		30mm	3-0	2	VIOLET	70		M3E30U70
			2-0	3	VIOLET	70		M2E30U70
			0	3.5	VIOLET	70		M0E30U70
			1	4	VIOLET	70		M1E30U70
E30U-D		30mm	3-0	2	VIOLET	70		M3E30U70D
			2-0	3	VIOLET	70		M2E30U70D




1/2 Circle		X-Trocar Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E37U  37mm		3-0	2	VIOLET	70		M3E37U70
		3-0	2	VIOLET	90		M3E37U90
		2-0	3	VIOLET	70		M2E37U70
		2-0	3	VIOLET	90		M2E37U90
		0	3.5	VIOLET	70		M0E37U70
		0	3.5	VIOLET	90		M0E37U90
		1	4	VIOLET	70		M1E37U70
		1	4	VIOLET	90		M1E37U90
E40U  40mm		0	3.5	VIOLET	90		M0E40U90
		1	4	VIOLET	90		M1E40U90
E48U  48mm		0	3.5	VIOLET	90		M0E48U90
		1	4	VIOLET	90		M1E48U90
		1	4	VIOLET	150		M1E48U150

3/8 Circle		X-Trocar Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E18C	 18mm	5-0	1	VIOLET	70		M5E18C70
		4-0	1.5	VIOLET	70		M4E18C70
E22C	 22mm	4-0	1.5	VIOLET	70		M4E22C70
		3-0	2	VIOLET	70		M3E22C70
E26C	 26mm	3-0	2	VIOLET	70		M3E26C70
		2-0	3	VIOLET	70		M2E26C70
E32C	 32mm	2-0	3	VIOLET	70		M2E32C70
		0	3.5	VIOLET	70		M0E32C70
		1	4	VIOLET	70		M1E32C70

1/2 , 5/8 Circle		Taper Point Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R26U-0-D	 26mm	3-0	2	UNDYED	70		M3R26U-070DU
		2-0	3	UNDYED	70		M2R26U-070DU

1/2 , 5/8 Circle		Taper Point Needles			◎		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R26U-0-D	 26mm	2-0	3	VIOLET	90		M2R26U-090D

1/2 Circle		Taper Point, Reverse Cutting Needles			◎ ▼		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R-C37U-D	 37mm	3-0	2	UNDYED	140		M3R-C37U140DU
		2-0	3	UNDYED	150	LOOP	M2R-C37U150DLU

1/2 Circle		Reverse Cutting, Taper Point Needles			▼ ◎		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C-R17U-D	 17mm	3-0	2	VIOLET	70		M3C-R17U70
C-R22U-D	 22mm	3-0	2	UNDYED	90		M3C-R22U90DU
C-R36U-D	 36mm	3-0	2	UNDYED	140		M3C-R36U140DU

Precut						
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
	5-0	1	VIOLET	3 X 45		M5TY345
	4-0	1.5	VIOLET	3 X 45		M4TY345
	4-0	1.5	VIOLET	5 X 70		M4TY570
	4-0	1.5	VIOLET	6 X 45		M4TY645
	3-0	2	VIOLET	1 X 140		M3TY1140
	3-0	2	VIOLET	3 X 45		M3TY345
	3-0	2	VIOLET	5 X 70		M3TY570
	3-0	2	VIOLET	6 X 45		M3TY645
	2-0	3	VIOLET	1 X 140		M2TY1140
	2-0	3	VIOLET	3 X 45		M2TY345
	2-0	3	VIOLET	5 X 70		M2TY570
	2-0	3	VIOLET	6 X 45		M2TY645
	2-0	3	UNDYED	6 X 45		M2TY645U
	0	3.5	VIOLET	1 X 140		M0TY1140
	0	3.5	VIOLET	3 X 45		M0TY345
	0	3.5	VIOLET	5 X 70		M0TY570

Precut							
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
	0	3.5	VIOLET	6 X 45		MOTY645	
	0	3.5	UNDYED	6 X 45		MOTY645U	
	1	4	VIOLET	1 X 140		M1TY1140	
	1	4	VIOLET	3 X 45		M1TY345	
	1	4	VIOLET	5 X 70		M1TY570	
	1	4	VIOLET	6 X 45		M1TY645	
	1	4	UNDYED	6 X 45		M1TY645U	