

Suture

Description: Monofilament For smooth passage through tissue.

Synthetic For minimal tissue reaction.
Absorbable By hydrolysis. Essentially complete in 182-238 days.

Tensile strength: Post implantation approximate % of original strength remaining

14 days	80% (3-0 and larger) 60% (4-0 and smaller)
28 days	70% (3-0 and larger) 40% (4-0 and smaller)
42 days	60% (3-0 and larger)
28 days	35% (4-0 and smaller)

Color: Violet or Undyed.



Range: 6-0 to 1 (USP). Supplied as needled sutures and ligatures.





Indications: PDO* sutures are intended for use in general soft tissue approximation, including in pediatric cardiovascular tissue, and in ophthalmic surgery (other than contact with cornea and sclera). PDO* sutures are particularly useful where the combination of absorbable suture and extended wound support is desirable. Typical areas of use include:



- Abdominal Wall Closure
- Intestinal Anastomoses
- Ligament and Tendon Repair


Full details in the instruction for use included in every package.












1/2 Circle		Taper Point Plus Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Q6.5-180U-150-D 6.5mm	 Length: 6.5mm Radius: 2.07mm Curvature: 180° Chord: 4.14mm Diameter: 150µ	7-0	0.5	VIOLET	70		D7Q6.5U70D
Q10-180U-330-D 10mm	 Length: 10mm Radius: 4.25mm Curvature: 180° Chord: 7.85mm Diameter: 330µ	6-0	0.7	VIOLET	70		D6Q10U70D


1/2 Circle		Taper Point Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R9.3-180U-200-D 9.3mm	 Length: 9.3mm Radius: 2.96mm Curvature: 180° Chord: 5.92mm Diameter: 200µ	6-0	0.7	VIOLET	70		D6R9.3U70D
R10-180U-330 10mm	 Length: 10mm Radius: 3.18mm Curvature: 180° Chord: 6.36mm Diameter: 330µ	6-0	0.7	VIOLET	75		D6R10U75
R11-180U-250 11mm	 Length: 11mm Radius: 3.5mm Curvature: 180° Chord: 7.00mm Diameter: 250µ	5-0	1	VIOLET	70		D5R11U70
R13-180U-350 13mm	 Length: 13mm Radius: 4.14mm Curvature: 180° Chord: 8.28mm Diameter: 350µ	6-0	0.7	UNDYED	45		D6R13U45U
		6-0	0.7	VIOLET	60		D6R13U60
		6-0	0.7	VIOLET	70		D6R13U70




1/2 Circle		Taper Point Needles			Micro Needles			
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R13-180U-350 13mm		Length: 13mm	Radius: 4.14mm	6-0	0.7	VIOLET	75	D6R13U75
		Curvature: 180°	Chord: 8.28mm	5-0	1	VIOLET	45	D5R13U45
		Diameter: 350µ		5-0	1	UNDYED	45	D5R13U45U
				5-0	1	UNDYED	60	D5R13U60
				5-0	1	VIOLET	70	D5R13U70
				5-0	1	VIOLET	75	D5R13U75
				5-0	1	VIOLET	90	D5R13U90
				4-0	1.5	VIOLET	45	D4R13U45
				4-0	1.5	VIOLET	70	D4R13U70
				4-0	1.5	VIOLET	75	D4R13U75
				3-0	2	VIOLET	70	HEAVY NEEDLE D3R13U70H
				2-0	3	VIOLET	70	D2R13U70
				2-0	3	VIOLET	70	HEAVY NEEDLE D2R13U70H
R13-180U-350-D 13mm		Length: 13mm	Radius: 4.14mm	6-0	0.7	VIOLET	70	D6R13U70D
		Curvature: 180°	Chord: 8.28mm	6-0	0.7	VIOLET	75	D6R13U75D
		Diameter: 350µ		5-0	1	VIOLET	70	D5R13U70D


1/2 Circle		Taper Point Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R13-180U-400-D 13mm		Length: 13mm	Radius: 4.14mm	5-0	1	VIOLET	70	D5R13U70D
		Curvature: 180°	Chord: 8.28mm					
		Diameter: 400µ						
		5-0	1	VIOLET	75		D5R13U75D	
		5-0	1	UNDYED	90		D5R13U90DU	



3/8 Circle		Taper Point Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R6.5-135C-130-D 6.5mm		Length: 6.5mm	Radius: 2.76mm	7-0	0.5	VIOLET	75	D7R6.5C75D
		Curvature: 135°	Chord: 5.1mm					
		Diameter: 130µ						
R8-135C-175-D 8mm		Length: 8mm	Radius: 3.4mm	7-0	0.5	VIOLET	60	D7R8C60D
		Curvature: 135°	Chord: 6.28mm					
		Diameter: 175µ						
R9-135C-200-D 9mm		Length: 9mm	Radius: 3.82mm	7-0	0.5	VIOLET	45	D7R9C45D
		Curvature: 135°	Chord: 7.06mm	7-0	0.5	VIOLET	60	D7R9C60D
			Diameter: 200µ	7-0	0.5	VIOLET	75	D7R9C75D
				6-0	0.7	VIOLET	60	D6R9C60D
				6-0	0.7	VIOLET	75	D6R9C75D



3/8 Circle		Taper Point Needles				Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R9.3-150C-175-B 9.3mm	 Length: 9.3mm Curvature: 150° Radius: 3.55mm Chord: 6.86mm Diameter: 175µ	7-0	0.5	VIOLET	70	X-BLACK*	D7R9.3C70B	
R9.3-135C-200-D 9.3mm	 Length: 9.3mm Curvature: 135° Radius: 3.95mm Chord: 7.3mm Diameter: 200µ	7-0	0.5	VIOLET	60		D7R9.3C60D	
		7-0	0.5	VIOLET	70		D7R9.3C70D	
		6-0	0.7	VIOLET	60		D6R9.3C60D	
		6-0	0.7	VIOLET	75		D6R9.3C75D	
		5-0	1	VIOLET	75		D5R9.3C75D	
R11-135C-330-D 11mm	 Length: 11mm Curvature: 135° Radius: 4.67mm Chord: 8.63mm Diameter: 330µ	6-0	0.7	VIOLET	45		D6R11C45D	
		6-0	0.7	VIOLET	70		D6R11C70D	
R13-135C-530 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 530µ	6-0	0.7	VIOLET	70		D6R13C70	
		5-0	1	VIOLET	70		D5R13C70	
R13-135C-530-D 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 530µ	7-0	0.5	VIOLET	45		D7R13C45D	
		6-0	0.7	VIOLET	45		D6R13C45D	
		6-0	0.7	VIOLET	60		D6R13C60D	
		6-0	0.7	VIOLET	70		D6R13C70D	
		5-0	1	VIOLET	70		D5R13C70D	



3/8 Circle		Taper Point Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R13-135C-530-D 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 530µ	5-0	1	VIOLET	90		D5R13C90D


3/8 Circle		Taper Cutting Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T9-135C-200 9mm	 Length: 9mm Curvature: 135° Radius: 3.82mm Chord: 7.06mm Diameter: 200µ	6-0	0.7	VIOLET	75		D6T9C75
T13-135C-350 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 350µ	6-0	0.7	VIOLET	70		D6T13C70
		5-0	1	VIOLET	70		D5T13C70
T13-135C-250-D 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 250µ	6-0	0.7	VIOLET	75		D6T13C75D
		5-0	1	VIOLET	75		D5T13C75D


3/8 Circle		X-CC Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
H9.3-135C-200 9.3mm	 Length: 9.3mm Curvature: 135° Radius: 3.95mm Chord: 7.3mm Diameter: 200µ	7-0	0.5	VIOLET	70		D7H9.3C70



3/8 Circle		X-CC Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
H9.3-135C-200-D 9.3mm		Length: 9.3mm	Radius: 3.95mm	VIOLET	70		D7H9.3C70D
		Curvature: 135°	Chord: 7.3mm				
H13-135C-300-D 13mm		Length: 13mm	Radius: 5.52mm	VIOLET	70		D6H13C70D
				Curvature: 135°	Chord: 10.2mm	6-0	0.7
		Diameter: 200µ	6-0	0.7	VIOLET	70	
		Diameter: 300µ		5-0	1		D5H13C70D




3/8 Circle		Reverse Cutting Needles				Micro Needles				
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code			
C12-135C-300 12mm		Length: 12mm	Radius: 5.1mm	VIOLET	45		D6C12C45			
				Curvature: 135°	Chord: 9.42mm	VIOLET	45		D5C12C45	
				Diameter: 300µ	5-0	1	VIOLET	45		D5C12C45
					5-0	1	VIOLET	60		D5C12C60
					5-0	1	VIOLET	75		D5C12C75
C13-135C-350 13mm		Length: 13mm	Radius: 5.52mm	VIOLET	45		D4C12C45			
				Curvature: 135°	Chord: 10.2mm	4-0	1.5	VIOLET	45	
		Diameter: 350µ		5-0	1		D5C13C75			



1/2 Circle		Reverse Cutting Precision Point Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC8-180U-200 8mm	 Length: 8mm Radius: 2.55mm Curvature: 180° Chord: 5.1mm Diameter: 200µ	5-0	1	VIOLET	45		D5PC8U45	
		5-0	1	UNDYED	45		D5PC8U45U	
		4-0	1.5	UNDYED	45		D4PC8U45U	
PC10-180U-250 10mm	 Length: 10mm Radius: 3.18mm Curvature: 180° Chord: 6.36mm Diameter: 250µ	6-0	0.7	UNDYED	45		D6PC10U45U	



3/8 Circle		Reverse Cutting Precision Point Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC12-135C-300 12mm	 Length: 12mm Radius: 5.1mm Curvature: 135° Chord: 9.42mm Diameter: 300µ	6-0	0.7	VIOLET	45		D6PC12C45	
		5-0	1	VIOLET	60		D5PC12C60	
		5-0	1	VIOLET	75		D5PC12C75	



1/2 Circle		X-Prime Reverse Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X8-180U-200 8mm	 Length: 8mm Radius: 2.55mm Curvature: 180° Chord: 5.1mm Diameter: 200µ	6-0	0.7	VIOLET	45		D6X8U45	



1/2 Circle		X-Prime Reverse Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X13-180U-350 13mm		Length: 13mm	Radius: 4.14mm	6-0	0.7	UNDYED	45	D6X13U45U
		Curvature: 180°	Chord: 8.28mm	4-0	1.5	VIOLET	45	D4X13U45
X13-180U-350-D 13mm		Length: 13mm	Radius: 4.14mm	6-0	0.7	VIOLET	75	D6X13U75D
		Curvature: 180°	Chord: 8.28mm	5-0	1	VIOLET	75	D5X13U75D



3/8 Circle		X-Prime Reverse Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X11-135C-250 11mm		Length: 11mm	Radius: 4.68mm	6-0	0.7	VIOLET	45	D6X11C45
		Curvature: 135°	Chord: 8.63mm	6-0	0.7	UNDYED	45	D6X11C45U
			Diameter: 250µ	5-0	1	UNDYED	45	D5X11C45U
X11-135C-250-D 11mm		Length: 11mm	Radius: 4.68mm	6-0	0.7	VIOLET	75	D6X11C75D
X13-135C-350 13mm		Length: 13mm	Radius: 5.52mm	5-0	1	VIOLET	45	D5X13C45
		Curvature: 135°	Chord: 10.2mm	5-0	1	UNDYED	45	D5X13C45U
			Diameter: 350µ	4-0	1.5	VIOLET	45	D4X13C45
				4-0	1.5	UNDYED	45	D4X13C45U



1/4 Circle		X-Prime Reverse Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X8-90Q-200-D		Length: 8mm Radius: 5.1mm Curvature: 90° Chord: 7.21mm Diameter: 200µ	6-0	0.7	VIOLET	45	D6X8Q45D
8mm							



3/8 Circle		Conventional Cutting Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K12-135C-300		Length: 12mm Radius: 5.1mm Curvature: 135° Chord: 9.42mm Diameter: 300µ	6-0	0.7	VIOLET	45	D6K12C45
12mm							




Straight		Conventional Cutting Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K13-0Y-380		Length: 13mm Radius: 0mm Curvature: 0° Chord: 13.00mm Diameter: 380µ	4-0	1.5	VIOLET	35	D4K13Y35
13mm							






3/8 Circle		Conventional Cutting Precision Point Needles 				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PK12-135C-300 12mm	 Length: 12mm Curvature: 135° Radius: 5.1mm Chord: 9.42mm Diameter: 300µ	6-0	0.7	VIOLET	45		D6PK12C45




3/8 Circle		X-Prime Cutting Needles 				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D13-135C-330 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 330µ	6-0	0.7	UNDYED	45		D6D13C45U





1/2 Circle		Conventional Spatula Needles 				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Z7.6-163U-200-D 7.6mm	 Length: 7.6mm Curvature: 163° Radius: 2.67mm Chord: 5.28mm Diameter: 200µ	7-0	0.5	VIOLET	45		D7Z7.6U45D




1/4 Circle		Conventional Spatula Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Z8-112Q-330-D	 8mm	Length: 8mm Curvature: 112°	Radius: 4.09mm Chord: 6.78mm Diameter: 330µ	6-0	0.7	UNDYED	45	D6Z8Q45DU


1/2 Circle		Micropoint Spatula Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
F5.5-160U-200	 5.5mm	Length: 5.5mm Curvature: 160°	Radius: 1.97mm Chord: 3.88mm Diameter: 200µ	7-0	0.5	VIOLET	45	D7F5.5U45
F5.5-160U-200-D	 5.5mm	Length: 5.5mm Curvature: 160°	Radius: 1.97mm Chord: 3.88mm Diameter: 200µ	7-0	0.5	VIOLET	20	D7F5.5U20D




1/2 Circle		Taper Point Plus Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q16U-H		16mm	3-0	2	VIOLET	70	HEAVY NEEDLE	D3Q16U70H
			2-0	3	VIOLET	70	HEAVY NEEDLE	D2Q16U70H
			0	3.5	VIOLET	70	HEAVY NEEDLE	D0Q16U70H
Q17U		17mm	5-0	1	VIOLET	70		D5Q17U70
			4-0	1.5	VIOLET	70		D4Q17U70
			3-0	2	VIOLET	70		D3Q17U70
Q17U-D		17mm	5-0	1	VIOLET	70	X-BLACK*	D5Q17U70DB
			5-0	1	VIOLET	70		D5Q17U70D
			4-0	1.5	VIOLET	70		D4Q17U70D
			3-0	2	VIOLET	70		D3Q17U70D
Q22U		22mm	4-0	1.5	VIOLET	70		D4Q22U70
			3-0	2	VIOLET	70		D3Q22U70
			2-0	3	VIOLET	70	HEAVY NEEDLE	D2Q22U70H
			0	3.5	VIOLET	70	HEAVY NEEDLE	D0Q22U70H
Q22U-D		22mm	4-0	1.5	VIOLET	70		D4Q22U70D
			4-0	1.5	VIOLET	70	X-BLACK*	D4Q22U70DB




1/2 Circle		Taper Point Plus Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Q22U-D-B		4-0	1.5	VIOLET	75	X-BLACK*	D4Q22U75DB
		4-0	1.5	VIOLET	4 X 70	CONTROL RELEASE*	D4Q22U70DB
		3-0	2	VIOLET	70		D3Q22U70D
		3-0	2	VIOLET	70	X-BLACK*	D3Q22U70DB
Q26U		4-0	1.5	VIOLET	70		D4Q26U70
		3-0	2	VIOLET	70		D3Q26U70
		3-0	2	VIOLET	70	X-BLACK*	D3Q26U70B
		2-0	3	VIOLET	70		D2Q26U70
		2-0	3	VIOLET	70	X-BLACK*	D2Q26U70B
		2-0	3	VIOLET	70	HEAVY NEEDLE	D2Q26U70H
		0	3.5	VIOLET	70	HEAVY NEEDLE	D0Q26U70H
		1	4	VIOLET	70	HEAVY NEEDLE	D1Q26U70H
Q26U-D		4-0	1.5	VIOLET	90	X-BLACK*	D4Q26U90DB
		3-0	2	VIOLET	70		D3Q26U70D
		3-0	2	VIOLET	90	X-BLACK*	D3Q26U90DB
		2-0	3	VIOLET	90		D2Q26U90D




1/2 Circle		Taper Point Plus Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q31U		31mm	3-0	2	VIOLET	70		D3Q31U70
			2-0	3	VIOLET	150	LOOP	D2Q31U150L
Q36U		36mm	3-0	2	VIOLET	70		D3Q36U70
			3-0	2	VIOLET	70	HEAVY NEEDLE	D3Q36U70H
			2-0	3	VIOLET	70		D2Q36U70
			2-0	3	VIOLET	70	HEAVY NEEDLE	D2Q36U70H
			2-0	3	VIOLET	150	HEAVY NEEDLE	D2Q36U150H
			0	3.5	VIOLET	70	HEAVY NEEDLE	D0Q36U70H
			0	3.5	VIOLET	90	HEAVY NEEDLE	D0Q36U90H
			1	4	VIOLET	70	HEAVY NEEDLE	D1Q36U70H
Q40U-H		40mm	2-0	3	VIOLET	90	HEAVY NEEDLE	D2Q40U90H
			0	3.5	VIOLET	90	HEAVY NEEDLE	D0Q40U90H
			0	3.5	VIOLET	150	HEAVY NEEDLE/LOOP	D0Q40U150HL
			1	4	VIOLET	90	HEAVY NEEDLE	D1Q40U90H
			1	4	VIOLET	150	HEAVY NEEDLE/LOOP	D1Q40U150HL



1/2 Circle		Taper Point Plus Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Q48U-H  48mm		0	3.5	VIOLET	90	HEAVY NEEDLE	D0Q48U90H
		0	3.5	VIOLET	150	HEAVY NEEDLE	D0Q48U150H
		0	3.5	VIOLET	150	HEAVY NEEDLE/LOOP	D0Q48U150HL
		0	3.5	VIOLET	150	HEAVY NEEDLE/LOOP	D0Q48U150HL
		1	4	VIOLET	90	HEAVY NEEDLE	D1Q48U90H
		1	4	VIOLET	150	HEAVY NEEDLE/LOOP	D1Q48U150HL
		1	4	VIOLET	240	HEAVY NEEDLE/LOOP	D1Q48U240HL
Q65U-H  65mm		0	3.5	VIOLET	150	HEAVY NEEDLE/LOOP	D0Q65U150HL
		1	4	VIOLET	150	HEAVY NEEDLE/LOOP	D1Q65U150HL
Q70U-H  70mm		0	3.5	VIOLET	70	HEAVY NEEDLE	D0Q70U70H




1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R16U  16mm		4-0	1.5	UNDYED	75		D4R16U75U


1/2 Circle		Taper Point Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R17U		17mm	5-0	1	VIOLET	60		D5R17U60
			5-0	1	VIOLET	75		D5R17U75
			4-0	1.5	VIOLET	45		D4R17U45
			4-0	1.5	VIOLET	60		D4R17U60
			4-0	1.5	UNDYED	75		D4R17U75U
			3-0	2	VIOLET	75		D3R17U75
R17U-D		17mm	5-0	1	VIOLET	75		D5R17U75D
			5-0	1	VIOLET	90		D5R17U90D
			4-0	1.5	VIOLET	70		D4R17U70D
			4-0	1.5	UNDYED	70		D4R17U70DU
			4-0	1.5	VIOLET	90		D4R17U90D
			3-0	2	VIOLET	90		D3R17U90D
			2-0	3	VIOLET	90		D2R17U90D
R20U		20mm	5-0	1	UNDYED	75		D5R20U75U
			4-0	1.5	VIOLET	75		D4R20U75
			3-0	2	VIOLET	75		D3R20U75


1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R20U-D		20mm	4-0	1.5	VIOLET	90		D4R20U90D
			3-0	2	VIOLET	90		D3R20U90D
R22U		22mm	4-0	1.5	VIOLET	60		D4R22U60
			4-0	1.5	VIOLET	70		D4R22U70
			3-0	2	VIOLET	45		D3R22U45
			3-0	2	VIOLET	70		D3R22U70
			3-0	2	VIOLET	120		D3R22U120
			0	3.5	VIOLET	90		D0R22U90
R22U-D		22mm	4-0	1.5	VIOLET	90		D4R22U90D
R26U		26mm	4-0	1.5	VIOLET	45		D4R26U45
			4-0	1.5	VIOLET	75		D4R26U75
			3-0	2	VIOLET	45		D3R26U45
			3-0	2	UNDYED	70		D3R26U70U
			3-0	2	VIOLET	75		D3R26U75
			3-0	2	VIOLET	90		D3R26U90



1/2 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R26U 	26mm	2-0	3	VIOLET	45	D2R26U45	
		2-0	3	VIOLET	75	D2R26U75	
		2-0	3	VIOLET	90	D2R26U90	
		2-0	3	VIOLET	120	D2R26U120	
		0	3.5	VIOLET	45	D0R26U45	
		0	3.5	VIOLET	120	D0R26U120	
R26U-D 	26mm	4-0	1.5	VIOLET	75	D4R26U75D	
		4-0	1.5	VIOLET	90	D4R26U90D	
		4-0	1.5	VIOLET	105	D4R26U105D	
		3-0	2	VIOLET	90	D3R26U90D	
		2-0	3	VIOLET	90	D2R26U90D	
R27U 	27mm	3-0	2	VIOLET	75	D3R27U75	
		2-0	3	VIOLET	75	D2R27U75	
		0	3.5	VIOLET	75	D0R27U75	
		0	3.5	VIOLET	75	HEAVY NEEDLE D0R27U75H	
		1	4	VIOLET	75	D1R27U75	



1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R30U		4-0	1.5	VIOLET	75		D4R30U75
		3-0	2	VIOLET	75		D3R30U75
		2-0	3	VIOLET	75		D2R30U75
		0	3.5	VIOLET	75		D0R30U75
	30mm						
R31U		3-0	2	VIOLET	75		D3R31U75
		2-0	3	VIOLET	75		D2R31U75
		0	3.5	VIOLET	75		D0R31U75
		1	4	VIOLET	75		D1R31U75
	31mm						
R36U		3-0	2	VIOLET	45		D3R36U45
		3-0	2	VIOLET	70		D3R36U70
		2-0	3	VIOLET	45		D2R36U45
		2-0	3	VIOLET	60		D2R36U60
		2-0	3	VIOLET	70		D2R36U70
		2-0	3	UNDYED	70		D2R36U70U
		2-0	3	VIOLET	75		D2R36U75
		2-0	3	VIOLET	90		D2R36U90
	36mm						





1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R36U 	36mm	2-0	3	VIOLET	150		D2R36U150
		0	3.5	VIOLET	45		D0R36U45
		0	3.5	VIOLET	60		D0R36U60
		0	3.5	VIOLET	70		D0R36U70
		0	3.5	VIOLET	90		D0R36U90
		0	3.5	VIOLET	150	LOOP	D0R36U150L
		1	4	VIOLET	45		D1R36U45
		1	4	VIOLET	60		D1R36U60
		1	4	VIOLET	70		D1R36U70
		1	4	VIOLET	90		D1R36U90
1	4	VIOLET	250		D1R36U250		
R36U-D 	36mm	2-0	3	VIOLET	90		D2R36U90D
R37U 	37mm	3-0	2	VIOLET	75		D3R37U75
		2-0	3	VIOLET	75		D2R37U75
		2-0	3	VIOLET	90		D2R37U90


1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R37U								
			0	3.5	VIOLET	75		D0R37U75
			0	3.5	VIOLET	90		D0R37U90
			0	3.5	VIOLET	75	HEAVY NEEDLE	D0R37U75H
			0	3.5	VIOLET	90	HEAVY NEEDLE	D0R37U90H
		37mm	1	4	VIOLET	75		D1R37U75
			1	4	VIOLET	75	HEAVY NEEDLE	D1R37U75H
			1	4	VIOLET	90		D1R37U90
			1	4	VIOLET	90	HEAVY NEEDLE	D1R37U90H
			1	4	VIOLET	150	LOOP	D1R37U150L
			2	5	VIOLET	75		DXR37U75
	2	5	VIOLET	90	HEAVY NEEDLE	DXR37U90H		
R40U								
			2-0	3	VIOLET	75		D2R40U75
			2-0	3	VIOLET	90		D2R40U90
			0	3.5	VIOLET	45		D0R40U45
		40mm	0	3.5	VIOLET	60		D0R40U60
	0	3.5	VIOLET	70		D0R40U70		






1/2 Circle		Taper Point Needles				
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R40U  40mm	0	3.5	VIOLET	75		D0R40U75
	0	3.5	VIOLET	90		D0R40U90
	0	3.5	VIOLET	90	HEAVY NEEDLE	D0R40U90H
	0	3.5	VIOLET	150		D0R40U150
	0	3.5	VIOLET	150	LOOP	D0R40U150L
	1	4	VIOLET	45		D1R40U45
	1	4	VIOLET	60		D1R40U60
	1	4	VIOLET	75		D1R40U75
	1	4	VIOLET	75	HEAVY NEEDLE	D1R40U75H
	1	4	VIOLET	90		D1R40U90
	1	4	VIOLET	90	HEAVY NEEDLE	D1R40U90H
	1	4	VIOLET	150	LOOP	D1R40U150L
	1	4	VIOLET	180		D1R40U180
	2	5	VIOLET	90		DXR40U90
	2	5	VIOLET	90	HEAVY NEEDLE	DXR40U90H


1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R48U 	48mm	3-0	2	VIOLET	90		D3R48U90
		2-0	3	VIOLET	90		D2R48U90
		2-0	3	VIOLET	150		D2R48U150
		0	3.5	VIOLET	45		D0R48U45
		0	3.5	VIOLET	60		D0R48U60
		0	3.5	VIOLET	90		D0R48U90
		0	3.5	VIOLET	150	LOOP	D0R48U150L
		0	3.5	VIOLET	180		D0R48U180
		1	4	VIOLET	45		D1R48U45
		1	4	VIOLET	60		D1R48U60
		1	4	VIOLET	70		D1R48U70
		1	4	VIOLET	90		D1R48U90
		1	4	VIOLET	150	LOOP	D1R48U150L
		2	5	VIOLET	90		DXR48U90
R48U-D 	48mm	0	3.5	VIOLET	45		D0R48U45D





1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R65U  65mm		0	3.5	VIOLET	150		D0R65U150
		0	3.5	VIOLET	250		D0R65U250
		1	4	VIOLET	75		D1R65U75
		1	4	VIOLET	100		D1R65U100
		1	4	VIOLET	150		D1R65U150
		1	4	VIOLET	180		D1R65U180
		1	4	VIOLET	250		D1R65U250
		2	5	VIOLET	90		DXR65U90
		2	5	VIOLET	150		DXR65U150
		2	5	VIOLET	250		DXR65U250
R70U  70mm		0	3.5	VIOLET	70		D0R70U70
		1	4	VIOLET	70		D1R70U70
		1	4	VIOLET	140		D1R70U140
		1	4	VIOLET	250		D1R70U250


3/8 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R17C		17mm	4-0	1.5	VIOLET	70		D4R17C70
			4-0	1.5	VIOLET	75		D4R17C75
R17C-D		17mm	5-0	1	VIOLET	75		D5R17C75D
			5-0	1	VIOLET	90		D5R17C90D
			4-0	1.5	VIOLET	75		D4R17C75D
R19C		19mm	3-0	2	VIOLET	70		D3R19C70
			3-0	2	UNDYED	70		D3R19C70U
R24C		24mm	2-0	3	UNDYED	70		D2R24C70U


5/8 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R260		26mm	3-0	2	VIOLET	70		D3R26070
			2-0	3	VIOLET	70		D2R26070
			0	3.5	VIOLET	70		D0R26070
			1	4	VIOLET	75		D1R26075


5/8 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R260-D 	26mm	3-0	2	VIOLET	70	D3R26070D	
R270 	27mm	2-0	3	VIOLET	75	D2R27075	
		0	3.5	VIOLET	45	D0R27045	
		0	3.5	VIOLET	75	D0R27075	
		1	4	VIOLET	75	D1R27075	
R270-D 	27mm	2-0	3	VIOLET	75	D2R27075D	
R360 	36mm	3-0	2	VIOLET	70	D3R36070	
		2-0	3	VIOLET	70	D2R36070	
R1100 	110mm	0	3.5	VIOLET	150	LOOP D0R1100150L	


Ski		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R19S		19mm	2-0	3	VIOLET	120	D2R19S120
			0	3.5	VIOLET	120	D0R19S120





Straight		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R19Y	 19mm	5-0	1	VIOLET	35		D5R19Y35
R19Y-D	 19mm	4-0	1.5	VIOLET	35		D4R19Y35D
R60Y	 60mm	3-0	2	VIOLET	75		D3R60Y75
		0	3.5	VIOLET	75		D0R60Y75
R65Y-D	 65mm	0	3.5	VIOLET	70		D0R65Y70D



1/2 Circle		X-GUARD Blunt Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
G45U	 45mm	1	4	VIOLET	90		D1G45U90


1/2 Circle		X-GUARD Blunt Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
G48U	 48mm	0	3.5	VIOLET	150	LOOP	DOG48U150L
		1	4	VIOLET	150	LOOP	D1G48U150L

5/8 Circle		X-GUARD Blunt Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
G1100	 110mm	0	3.5	VIOLET	240	LOOP	DOG1100240L



1/2 Circle		Blunt Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
B22U	 22mm	3-0	2	VIOLET	45		D3B22U45


1/2 Circle		Blunt Point Needles			○		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
B26U		3-0	2	VIOLET	45		D3B26U45
		3-0	2	VIOLET	70		D3B26U70
		2-0	3	VIOLET	45		D2B26U45
		2-0	3	VIOLET	70		D2B26U70
		0	3.5	VIOLET	70		D0B26U70
B36U		2-0	3	VIOLET	70		D2B36U70
		0	3.5	VIOLET	45		D0B36U45
		0	3.5	VIOLET	70		D0B36U70
		0	3.5	VIOLET	90		D0B36U90
		1	4	VIOLET	45		D1B36U45
		1	4	VIOLET	90		D1B36U90
B37U		2-0	3	VIOLET	90		D2B37U90
		0	3.5	VIOLET	90		D0B37U90
B40U		2-0	3	VIOLET	90		D2B40U90
		0	3.5	VIOLET	90		D0B40U90
		1	4	VIOLET	90		D1B40U90


1/2 Circle		Blunt Point Needles			○		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
B48U  48mm		0	3.5	VIOLET	45		D0B48U45
		0	3.5	VIOLET	90		D0B48U90
		0	3.5	VIOLET	150	LOOP	D0B48U150L
		1	4	VIOLET	45		D1B48U45
		1	4	VIOLET	90		D1B48U90
		1	4	VIOLET	150	LOOP	D1B48U150L
		1	4	VIOLET	250		D1B48U250
B65U  65mm		0	3.5	VIOLET	90		D0B65U90
		0	3.5	VIOLET	150		D0B65U150
		1	4	VIOLET	140		D1B65U140
		1	4	VIOLET	250		D1B65U250


1/2 Circle		Taper Cutting Needles			◊		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T17U  17mm		5-0	1	VIOLET	70		D5T17U70
		4-0	1.5	VIOLET	70		D4T17U70




1/2 Circle		Taper Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T17U	 17mm	3-0	2	VIOLET	70		D3T17U70
		2-0	3	VIOLET	70		D2T17U70
T17U-D	 17mm	4-0	1.5	VIOLET	90		D4T17U90D
T22U	 22mm	4-0	1.5	VIOLET	70		D4T22U70
		3-0	2	VIOLET	70		D3T22U70
T26U	 26mm	4-0	1.5	VIOLET	70		D4T26U70
		3-0	2	VIOLET	70		D3T26U70
		2-0	3	VIOLET	70		D2T26U70
T30U-D	 30mm	0	3.5	VIOLET	45		D0T30U45D
T36U	 36mm	2-0	3	VIOLET	70		D2T36U70
		0	3.5	VIOLET	90		D0T36U90
		1	4	VIOLET	90		D1T36U90
T37U	 37mm	0	3.5	VIOLET	90		D0T37U90




1/2 Circle		Taper Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
T40U		40mm	0	3.5	VIOLET	150	LOOP	DOT40U150L
			1	4	VIOLET	150	LOOP	D1T40U150L
T48U		48mm	0	3.5	VIOLET	150	LOOP	DOT48U150L
			1	4	VIOLET	90		D1T48U90
			1	4	VIOLET	150	LOOP	D1T48U150L





Hook		Taper Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
T31J		31mm	0	3.5	VIOLET	70		DOT31J70
			1	4	VIOLET	70		D1T31J70


Straight		Taper Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
T38Y-D		38mm	4-0	1.5	VIOLET	35		D4T38Y35D



5/8 Circle		X-CC Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
 <p>63mm</p>	1	4	VIOLET	150	LOOP	D1H630150L	


1/2 Circle		Reverse Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
 <p>20mm</p>	3-0	2	VIOLET	45		D3C20U45	
 <p>22mm</p>	4-0	1.5	VIOLET	45		D4C22U45	
	3-0	2	UNDYED	70	HEAVY NEEDLE	D3C22U70HU	
	1	4	VIOLET	70		D1C22U70	
	1	4	VIOLET	70	HEAVY NEEDLE	D1C22U70H	
 <p>26mm</p>	3-0	2	VIOLET	45		D3C26U45	
	3-0	2	VIOLET	75		D3C26U75	
	2-0	3	VIOLET	45		D2C26U45	
	2-0	3	UNDYED	70		D2C26U70U	








1/2 Circle		Reverse Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C26U 	26mm	2-0	3	VIOLET	70	HEAVY NEEDLE	D2C26U70H
		2-0	3	VIOLET	75		D2C26U75
		0	3.5	VIOLET	45		D0C26U45
		0	3.5	VIOLET	70		D0C26U70
		0	3.5	VIOLET	70	HEAVY NEEDLE	D0C26U70H
		1	4	VIOLET	45		D1C26U45
		1	4	VIOLET	70	HEAVY NEEDLE	D1C26U70H
		1	4	VIOLET	75		D1C26U75
C26U-D 	26mm	2	5	VIOLET	45		DXC26U45D
C36U 	36mm	2-0	3	VIOLET	70	HEAVY NEEDLE	D2C36U70H
		0	3.5	VIOLET	45		D0C36U45
		0	3.5	VIOLET	60		D0C36U60
		0	3.5	VIOLET	70	HEAVY NEEDLE	D0C36U70H
		1	4	VIOLET	45		D1C36U45
		1	4	VIOLET	60		D1C36U60



1/2 Circle		Reverse Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C36U-H		36mm	1	4	VIOLET	70	HEAVY NEEDLE	D1C36U70H
			C37U		37mm	2-0	3	VIOLET
0	3.5	VIOLET	75			D0C37U75		
1	4	VIOLET	75			D1C37U75		
C40U		40mm	0	3.5	VIOLET	60		D0C40U60
			0	3.5	VIOLET	90	HEAVY NEEDLE	D0C40U90H
			1	4	VIOLET	45		D1C40U45
			1	4	VIOLET	60		D1C40U60
			1	4	VIOLET	75		D1C40U75
			1	4	VIOLET	90	HEAVY NEEDLE	D1C40U90H
			1	4	VIOLET	250		D1C40U250
			2	5	VIOLET	70		DXC40U70
			2	5	VIOLET	70	HEAVY NEEDLE	DXC40U70H
			C48U		48mm	3-0	2	VIOLET
2-0	3	VIOLET				150		D2C48U150




1/2 Circle		Reverse Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C48U  48mm	0	3.5	VIOLET	45		D0C48U45	
	0	3.5	VIOLET	90	HEAVY NEEDLE	D0C48U90H	
	0	3.5	VIOLET	150		D0C48U150	
	1	4	VIOLET	70		D1C48U70	
	1	4	VIOLET	90		D1C48U90	
	1	4	VIOLET	90	HEAVY NEEDLE	D1C48U90H	




3/8 Circle		Reverse Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C16C  16mm	5-0	1	VIOLET	75		D5C16C75	
	4-0	1.5	VIOLET	75		D4C16C75	
C19C  19mm	5-0	1	VIOLET	45		D5C19C45	
	5-0	1	UNDYED	45		D5C19C45U	
	5-0	1	VIOLET	60		D5C19C60	
	5-0	1	UNDYED	70		D5C19C70U	
	5-0	1	UNDYED	75		D5C19C75U	



3/8 Circle		Reverse Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C19C 	19mm	4-0	1.5	VIOLET	30	LOOP	D4C19C30L
		4-0	1.5	VIOLET	45		D4C19C45
		4-0	1.5	UNDYED	45	SLIM	D4C19C45SU
		4-0	1.5	UNDYED	70		D4C19C70U
		4-0	1.5	UNDYED	70	SLIM	D4C19C70SU
		4-0	1.5	VIOLET	75		D4C19C75
		3-0	2	UNDYED	45		D3C19C45U
		3-0	2	UNDYED	70		D3C19C70U
		3-0	2	VIOLET	75		D3C19C75
C24C 	24mm	4-0	1.5	UNDYED	45		D4C24C45U
		4-0	1.5	UNDYED	70		D4C24C70U
		4-0	1.5	VIOLET	70		D4C24C70
		4-0	1.5	VIOLET	75		D4C24C75
		3-0	2	VIOLET	70		D3C24C70
		3-0	2	UNDYED	70		D3C24C70U
		3-0	2	VIOLET	75		D3C24C75

3/8 Circle		Reverse Cutting Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
C24C		24mm	2-0	3	UNDYED	70		D2C24C70U
			2-0	3	VIOLET	75		D2C24C75
C26C		26mm	3-0	2	VIOLET	75		D3C26C75
			2-0	3	VIOLET	45		D2C26C45
			2-0	3	VIOLET	75		D2C26C75
C30C		30mm	3-0	2	VIOLET	75		D3C30C75
			2-0	3	VIOLET	75		D2C30C75
C31C		31mm	2-0	3	VIOLET	75		D2C31C75
C36C		36mm	2-0	3	VIOLET	75		D2C36C75
C37C		37mm	0	3.5	VIOLET	75		D0C37C75
			1	4	VIOLET	75		D1C37C75
C40C		40mm	0	3.5	VIOLET	75		D0C40C75
			1	4	VIOLET	90		D1C40C90
			1	4	VIOLET	90	HEAVY NEEDLE	D1C40C90H





Straight		Reverse Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
C60Y 	60mm	4-0	1.5	VIOLET	70		D4C60Y70
		4-0	1.5	UNDYED	70		D4C60Y70U
		3-0	2	VIOLET	70		D3C60Y70
		3-0	2	UNDYED	70		D3C60Y70U
		2-0	3	VIOLET	70		D2C60Y70
C60Y-D 	60mm	2-0	3	VIOLET	45		D2C60Y45D
		1	4	VIOLET	90		D1C60Y90D


3/8 Circle		Reverse Cutting Precision Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PC16C 	16mm	5-0	1	VIOLET	75		D5PC16C75
PC19C 	19mm	5-0	1	VIOLET	60		D5PC19C60
		4-0	1.5	VIOLET	75		D4PC19C75
PC26C 	26mm	3-0	2	VIOLET	75		D3PC26C75
		2-0	3	VIOLET	75		D2PC26C75


3/8 Circle		Reverse Cutting Precision Point Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PC31C	 31mm	2-0	3	VIOLET	75		D2PC31C75
PC36C	 36mm	2-0	3	VIOLET	75		D2PC36C75
PC37C	 37mm	0	3.5	VIOLET	75		D0PC37C75


1/2 Circle		X-Prime Reverse Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X16U	 16mm	4-0	1.5	VIOLET	45		D4X16U45
		4-0	1.5	UNDYED	45		D4X16U45U
		3-0	2	UNDYED	45		D3X16U45U
X26U	 26mm	3-0	2	VIOLET	70		D3X26U70
		2-0	3	VIOLET	70		D2X26U70










3/8 Circle		X-Prime Reverse Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X16C		16mm	5-0	1	UNDYED	45		D5X16C45U
			4-0	1.5	UNDYED	45		D4X16C45U
X19C		19mm	5-0	1	UNDYED	45		D5X19C45U
			4-0	1.5	VIOLET	45		D4X19C45
			4-0	1.5	UNDYED	45		D4X19C45U
			4-0	1.5	VIOLET	45	SLIM	D4X19C45S
			4-0	1.5	UNDYED	45	SLIM	D4X19C45SU
			4-0	1.5	VIOLET	70		D4X19C70
			3-0	2	VIOLET	45		D3X19C45
			3-0	2	UNDYED	45		D3X19C45U
			3-0	2	VIOLET	70		D3X19C70
X24C		24mm	4-0	1.5	UNDYED	45		D4X24C45U
			4-0	1.5	VIOLET	60		D4X24C60
			3-0	2	UNDYED	45		D3X24C45U
			3-0	2	VIOLET	70		D3X24C70
			2-0	3	UNDYED	45		D2X24C45U






3/8 Circle		X-Prime Reverse Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X26C  26mm	4-0	1.5	VIOLET	45		D4X26C45	
	3-0	2	VIOLET	45		D3X26C45	
	3-0	2	VIOLET	70		D3X26C70	
	2-0	3	UNDYED	45		D2X26C45U	
	2-0	3	UNDYED	70		D2X26C70U	


1/2 Circle		Conventional Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
K60U  60mm	1	4	VIOLET	250		D1K60U250	
	2	5	VIOLET	70		DXK60U70	
	2	5	VIOLET	150		DXK60U150	





3/8 Circle		Conventional Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
K16C  16mm	5-0	1	VIOLET	75		D5K16C75	
	4-0	1.5	VIOLET	75		D4K16C75	



3/8 Circle		Conventional Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
K19C		19mm	5-0	1	VIOLET	75		D5K19C75
			4-0	1.5	VIOLET	75		D4K19C75
K26C		26mm	3-0	2	VIOLET	75		D3K26C75
			2-0	3	VIOLET	75		D2K26C75
K30C		30mm	3-0	2	VIOLET	75		D3K30C75
K31C		31mm	2-0	3	VIOLET	75		D2K31C75
K36C		36mm	2-0	3	VIOLET	75		D2K36C75

Straight		Conventional Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
K60Y		60mm	4-0	1.5	VIOLET	75		D4K60Y75
			3-0	2	VIOLET	75		D3K60Y75
K254Y-D		254mm	2-0	3	VIOLET	75		D2K254Y75D

3/8 Circle		Conventional Cutting Precision Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PK16C		16mm	5-0	1	VIOLET	75	D5PK16C75
			4-0	1.5	VIOLET	75	D4PK16C75
PK19C		19mm	5-0	1	VIOLET	75	D5PK19C75
			4-0	1.5	VIOLET	75	D4PK19C75
PK26C		26mm	3-0	2	VIOLET	75	D3PK26C75
			2-0	3	VIOLET	75	D2PK26C75
PK30C		30mm	3-0	2	VIOLET	75	D3PK30C75
PK31C		31mm	2-0	3	VIOLET	75	D2PK31C75

1/2 Circle		X-Prime Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D17U-D		17mm	5-0	1	VIOLET	75	D5D17U75D

3/8 Circle		X-Prime Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D16C 	16mm	6-0	0.7	VIOLET	45		D6D16C45
		5-0	1	UNDYED	45		D5D16C45U
		4-0	1.5	UNDYED	45		D4D16C45U
D19C 	19mm	4-0	1.5	UNDYED	45		D4D19C45U
		3-0	2	UNDYED	45		D3D19C45U
D36C 	36mm	2-0	3	UNDYED	70		D2D36C70U

Straight		X-Trocar Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E254Y-D 	254mm	2-0	3	VIOLET	70		D2E254Y70D

Precut						
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
	3-0	2	VIOLET	3 X 45		D3TY345
	2-0	3	VIOLET	3 X 45		D2TY345
	0	3.5	VIOLET	150		D0TY150
	0	3.5	VIOLET	3 X 45		D0TY345