

Suture Description:

Braided	For easy handling and secure knot tying.
Coated	For smooth passage through tissue and easy knot tie down.
Synthetic	For minimal tissue reaction.
Absorbable	By hydrolysis. Essentially complete in 60-90 days

Tensile strength:




Post implantation approximate % original strength remaining	
14 days	75%
21 days	50% (6-0 and larger) 40% (7-0 and smaller)
28 days	25% (6-0 and larger)




Color: Violet or Undyed.



Range: 10-0 to 6 (USP). Supplied as needled sutures and ligatures.
Sizes 10-0 and 9-0 are monofilament in structure

Indications: VIPRONE* sutures are intended for use in soft tissue approximation and /or ligation, including use in ophthalmic surgery, peripheral nerve adaptation and microsurgery for vessels less than 2 mm in diameter.
Typical areas of use include:
Fascia Closure – Due to the 28-Day Claim
Subcutaneous Fat
Joint Capsule
Uterus
Full details in the instruction for use included in every package.









1/2 Circle		Taper Point Plus Needles				Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q6.5-180U-150 6.5mm		Length: 6.5mm	Radius: 2.07mm	8-0	0.4	VIOLET	30	G8Q6.5U30
		Curvature: 180°	Chord: 4.14mm	7-0	0.5	VIOLET	30	G7Q6.5U30
R10-180U-300 10mm		Length: 10mm	Radius: 3.18mm	5-0	1	VIOLET	70	G5Q10U70
		Curvature: 180°	Chord: 6.36mm					
Q13-180U-350 13mm		Length: 13mm	Radius: 4.14mm	6-0	0.7	UNDYED	45	G6Q13U45U
				6-0	0.7	VIOLET	70	G6Q13U70
				4-0	1.5	UNDYED	45	G4Q13U45U

1/2 Circle		Taper Point Needles				Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R7-180U-130 7mm		Length: 7mm	Radius: 1.91mm	7-0	0.5	VIOLET	30	G7R6U30
				7-0	0.5	VIOLET	45	G7R6U45
R10-180U-250 10mm		Length: 10mm	Radius: 3.18mm	6-0	0.7	VIOLET	70	G6R10U70
				5-0	1	VIOLET	70	G5R10U70
R12-180U-400 12mm		Length: 12mm	Radius: 3.82mm	6-0	0.7	VIOLET	45	G6R12U45
				6-0	0.7	UNDYED	45	G6R12U45U
				5-0	1	UNDYED	45	G5R12U45U



1/2 Circle		Taper Point Needles			Micro Needles			
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R12-180U-400 12mm		Length: 12mm Curvature: 180°	Radius: 3.82mm Chord: 7.64mm Diameter: 400µ	4-0	1.5	UNDYED	45	G4R12U45U
				6-0	0.7	VIOLET	45	G6R13U45
R13-180U-400 13mm		Length: 13mm Curvature: 180°	Radius: 4.14mm Chord: 8.28mm Diameter: 400µ	6-0	0.7	UNDYED	45	G6R13U45U
				6-0	0.7	VIOLET	60	G6R13U60
				6-0	0.7	UNDYED	60	G6R13U60U
				6-0	0.7	VIOLET	70	G6R13U70
				6-0	0.7	VIOLET	75	G6R13U75
				5-0	1	VIOLET	45	G5R13U45
				5-0	1	UNDYED	45	G5R13U45U
				5-0	1	VIOLET	70	G5R13U70
				5-0	1	UNDYED	75	G5R13U75U
				4-0	1.5	VIOLET	45	G4R13U45
				4-0	1.5	VIOLET	70	G4R13U70
				4-0	1.5	VIOLET	75	G4R13U75
				4-0	1.5	UNDYED	75	G4R13U75U








1/2 Circle		Taper Point Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R13-180U-400	 <p>Length: 13mm Curvature: 180°</p> <p>Radius: 4.14mm Chord: 8.28mm Diameter: 400µ</p>	3-0	2	VIOLET	70		G3R13U70
13mm							



3/8 Circle		Taper Point Needles				Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R4.8-100C-75	 <p>Length: 4.8mm Curvature: 100°</p> <p>Radius: 2.75mm Chord: 4.21mm Diameter: 75µ</p>	10-0	0.2	VIOLET	10	MONOFILAMENT	GAR4.8C10
4.8mm							
R4.9-100C-75	 <p>Length: 4.9mm Curvature: 100°</p> <p>Radius: 2.81mm Chord: 4.31mm Diameter: 75µ</p>	10-0	0.2	VIOLET	15	MONOFILAMENT	GAR4.9C15
4.9mm							
R6-135C-140	 <p>Length: 6mm Curvature: 135°</p> <p>Radius: 2.55mm Chord: 4.71mm Diameter: 140µ</p>	8-0	0.4	VIOLET	13		G8R6C13
6mm		8-0	0.4	VIOLET	20	HEAVY NEEDLE	G8R6C20H
R6-135C-140-D	 <p>Length: 6mm Curvature: 135°</p> <p>Radius: 2.55mm Chord: 4.71mm Diameter: 140µ</p>	8-0	0.4	VIOLET	30		G8R6C30D
6mm							
R6.5-135C-130	 <p>Length: 6.5mm Curvature: 135°</p> <p>Radius: 2.76mm Chord: 5.1mm Diameter: 130µ</p>	9-0	0.3	VIOLET	15	MONOFILAMENT	G9R6.5C15
6.5mm		8-0	0.4	VIOLET	20		G8R6.5C20


3/8 Circle		Taper Point Needles					Micro Needles	
	Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R8-135C-230 8mm	 Length: 8mm Curvature: 135° Radius: 3.4mm Chord: 6.28mm Diameter: 230µ		8-0	0.4	VIOLET	30		G8R8C30
R9.3-135C-200 9.3mm	 Length: 9.3mm Curvature: 135° Radius: 3.95mm Chord: 7.3mm Diameter: 200µ		8-0	0.4	VIOLET	30		G8R9.3C30
			6-0	0.7	VIOLET	45		G6R9.3C45
R10-135C-200 10mm	 Length: 10mm Curvature: 135° Radius: 4.25mm Chord: 7.85mm Diameter: 200µ		7-0	0.5	VIOLET	45		G7R10C45
			6-0	0.7	VIOLET	45		G6R10C45
			5-0	1	VIOLET	45		G5R10C45
R12-135C-230 12mm	 Length: 12mm Curvature: 135° Radius: 5.1mm Chord: 9.42mm Diameter: 230µ		6-0	0.7	VIOLET	70		G6R12C70
			5-0	1	VIOLET	45		G5R12C45
			5-0	1	VIOLET	70		G5R12C70
			4-0	1.5	VIOLET	60		G4R12C60
			4-0	1.5	VIOLET	70		G4R12C70
R12-135C-230-D 12mm	 Length: 12mm Curvature: 135° Radius: 5.1mm Chord: 9.42mm Diameter: 230µ		6-0	0.7	UNDYED	45		G6R12C45DU
			5-0	1	UNDYED	45		G5R12C45DU
R13-135C-530 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 530µ		6-0	0.7	VIOLET	75		G6R13C75
			5-0	1	VIOLET	75		G5R13C75


1/2 Circle		Taper Point Plus (Heavier Wire Diameter) 					Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code		
R13-180U-480-H 13mm		Length: 13mm	Radius: 4.14mm	3-0	2	VIOLET	70	HEAVY NEEDLE	G3R13U70H
		Curvature: 180°	Chord: 8.28mm	0	3.5	VIOLET	70	HEAVY NEEDLE	G0R13U70H
		Diameter: 480µ							


3/8 Circle		Taper Cutting Needles 					Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code		
T12-135C-250 12mm		Length: 12mm	Radius: 5.1mm	6-0	0.7	VIOLET	70		G6T12C70
		Curvature: 135°	Chord: 9.42mm	5-0	1	VIOLET	70		G5T12C70
			Diameter: 250µ	4-0	1.5	VIOLET	70		G4T12C70
T13-135C-350 13mm		Length: 13mm	Radius: 5.52mm	6-0	0.7	VIOLET	70		G6T13C70
		Curvature: 135°	Chord: 10.2mm	5-0	1	VIOLET	70		G5T13C70
		Diameter: 350µ							


5/8 Circle		Taper Cutting Needles 					Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code		
T10-2250-330 10mm		Length: 10mm	Radius: 2.55mm	6-0	0.7	VIOLET	30		G6T10030
		Curvature: 225°	Chord: 4.71mm						
		Diameter: 330µ							



3/8 Circle		Reverse Cutting Precision Point Needles				Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC11-135C-350 11mm		Length: 11mm Curvature: 135°	Radius: 4.67mm Chord: 8.63mm Diameter: 350µ	5-0	1	UNDYED	45	G5PC11C45U
				6-0	0.7	UNDYED	75	G6PC12C75U
PC12-135C-300 12mm		Length: 12mm Curvature: 135°	Radius: 5.1mm Chord: 9.42mm Diameter: 300µ	6-0	0.7	UNDYED	45	G6PC13C45U
				5-0	1	VIOLET	45	G5PC13C45
				5-0	1	UNDYED	45	G5PC13C45U
				4-0	1.5	VIOLET	45	G4PC13C45
				4-0	1.5	UNDYED	45	G4PC13C45U





3/8 Circle		X-Prime Reverse Needles				Micro Needles		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X11-135C-250 11mm		Length: 11mm Curvature: 135°	Radius: 4.67mm Chord: 8.63mm Diameter: 250µ	6-0	0.7	UNDYED	45	G6X11C45U
				5-0	1	UNDYED	45	G5X11C45U



3/8 Circle			X-Prime Reverse Needles				Micro Needles	
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X13-135C-350 13mm		Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 350µ	6-0	0.7	UNDYED	45		G6X13C45U
			5-0	1	UNDYED	45		G5X13C45U
			4-0	1.5	UNDYED	45		G4X13C45U



3/8 Circle			Conventional Cutting Needles				Micro Needles	
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K11-135C-300 11mm		Length: 11mm Curvature: 135° Radius: 4.67mm Chord: 8.63mm Diameter: 300µ	6-0	0.7	UNDYED	45		G6K11C45U
			5-0	1	UNDYED	45		G5K11C45U
			4-0	1.5	UNDYED	45		G4K11C45U



Straight			Conventional Cutting Needles				Micro Needles	
Needle			Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K13-0Y-380-D 13mm		Length: 13mm Curvature: 0° Radius: 0mm Chord: 13mm Diameter: 380µ	4-0	1.5	VIOLET	45		G4K13Y45D




3/8 Circle		Conventional Cutting Precision Point Needles 					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PK13-135C-380 13mm	 Length: 13mm Curvature: 135° Radius: 5.52mm Chord: 10.2mm Diameter: 380µ	5-0	1	UNDYED	45		G5PK13C45U	




3/8 Circle		Reverse Cutting Micropoint Needles 					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
M6.2-135C-200-D 6.2mm	 Length: 6.2mm Curvature: 135° Radius: 2.63mm Chord: 4.86mm Diameter: 200µ	10-0	0.2	VIOLET	30	MONOFILAMENT	GAM6.2C30D	
M11-135C-250-D 11mm	 Length: 11mm Curvature: 135° Radius: 4.67mm Chord: 8.63mm Diameter: 250µ	6-0	0.7	VIOLET	45		G6M11C45D	
M13-155C-350-D 13mm	 Length: 13mm Curvature: 155° Radius: 4.81mm Chord: 9.39mm Diameter: 350µ	5-0	1	VIOLET	45		G5M13C45D	



1/4 Circle		Reverse Cutting Micropoint Needles 					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
M8-110Q-200-D 8mm	 Length: 8mm Curvature: 110° Radius: 4.17mm Chord: 6.83mm Diameter: 200µ	6-0	0.7	VIOLET	45		G6M8Q45D	

1/4 Circle		Center Point Spatula Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
08.5-90Q-280-D 8.5mm		Length: 8.5mm	Radius: 5.00mm	6-0	0.7	VIOLET	45	G608.5Q45D
		Curvature: 90°	Chord: 7.65mm	5-0	1	VIOLET	45	G508.5Q45D
		Diameter: 280µ						



1/2 Circle		Conventional Spatula Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Z8-170U-300-D 8mm		Length: 8mm	Radius: 2.7 mm	6-0	0.7	VIOLET	45	G6Z8U45D
		Curvature: 170°	Chord: 5.38mm	6-0	0.7	UNDYED	45	G6Z8U45DU
		Diameter: 300µ		5-0	1	VIOLET	45	G5Z8U45D






3/8 Circle		Conventional Spatula Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Z6.4-135C-200-D 6.4mm		Length: 6.4mm	Radius: 2.72 mm	7-0	0.5	UNDYED	30	G7Z6.4C30DU
		Curvature: 135°	Chord: 5.03mm	6-0	0.7	UNDYED	45	G6Z6.4C45DU
Z11-135C-240-D 11mm		Length: 11mm	Radius: 4.67 mm	5-0	1	UNDYED	45	G5Z11C45DU
		Curvature: 135°	Chord: 8.63mm					
		Diameter: 240µ						







1/4 Circle		Conventional Spatula Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Z8-112Q-330 8mm	 Length: 8mm Curvature: 112° Radius: 4.09mm Chord: 6.78mm Diameter: 330µ	6-0	0.7	UNDYED	45		G6Z8Q45DU	
		6-0	0.7	VIOLET	45	X-BLACK	G6Z8Q45DB	
Z8-112Q-330-D 8mm	 Length: 8mm Curvature: 112° Radius: 4.09mm Chord: 6.78mm Diameter: 330µ	6-0	0.7	VIOLET	20		G6Z8Q20D	
		6-0	0.7	VIOLET	45		G6Z8Q45D	
		6-0	0.7	VIOLET	45	SLIM	G6Z8Q45DS	
		5-0	1	VIOLET	45		G5Z8Q45D	
		5-0	1	UNDYED	45		G5Z8Q45DU	
		5-0	1	VIOLET	45	X-BLACK	G5Z8Q45DB	
		4-0	1.5	VIOLET	45		G4Z8Q45D	



1/2 Circle		Micropoint Spatula Needles					Micro Needles	
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
F5.5-180U-100-D 5.5mm	 Length: 5.5mm Curvature: 180° Radius: 1.75mm Chord: 3.5mm Diameter: 100µ	10-0	0.2	VIOLET	30		GAF5.5U30D	
		9-0	0.3	VIOLET	30		G9F5.5U30D	
		8-0	0.4	VIOLET	20		G8F5.5U20D	







1/2 Circle		Micropoint Spatula Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
F7-175U-200-D 7mm	 Length: 7mm Curvature: 175°	Radius: 2.29mm Chord: 4.58mm Diameter: 200µ	7-0	0.5	VIOLET	45		G7F7U45D




3/8 Circle		Micropoint Spatula Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
F6-135C-140-D 6mm	 Length: 6mm Curvature: 135°	Radius: 2.55mm Chord: 4.71mm Diameter: 140µ	8-0	0.4	VIOLET	30		G8F6C30D
F6.5-140C-150 6.5mm	 Length: 6.5mm Curvature: 140°	Radius: 2.66mm Chord: 5.00mm Diameter: 150µ	8-0	0.4	VIOLET	30		G8F6.5C30
			7-0	0.5	VIOLET	30		G7F6.5C30
F6.5-140C-150-D 6.5mm	 Length: 6.5mm Curvature: 140°	Radius: 2.66mm Chord: 5.00mm Diameter: 150µ	10-0	0.2	VIOLET	30	MONOFILAMENT	GAF6.5C30D
			9-0	0.3	VIOLET	30		G9F6.5C30D
			9-0	0.3	VIOLET	30	MONOFILAMENT	G9F6.5C30D
			8-0	0.4	VIOLET	20		G8F6.5C20D
			8-0	0.4	VIOLET	30		G8F6.5C30D
F7-0	0.5	VIOLET	30		G7F6.5C30D			
F8-135C-200-D 8mm	 Length: 8mm Curvature: 135°	Radius: 3.04mm Chord: 6.28mm Diameter: 200µ	7-0	0.5	VIOLET	45		G7F8C45D



1/4 Circle		Micropoint Spatula Needles						Micro Needles
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
F6.5-90Q-240-D 6.5mm	 Length: 6.5mm Radius: 4.14mm Curvature: 90° Chord: 5.85mm Diameter: 240µ	6-0	0.7	VIOLET	45		G6F6.5Q45D	
F6.6-90Q-200-D 6.6mm	 Length: 6.6mm Radius: 4.14mm Curvature: 90° Chord: 5.85mm Diameter: 200µ	6-0	0.7	VIOLET	45		G6F6.6Q45D	
F7.62-90Q-200-D 7.62mm	 Length: 7.62mm Radius: 4.84mm Curvature: 90° Chord: 6.84mm Diameter: 200µ	6-0	0.5	VIOLET	30		G6F7.62Q30D	
		6-0	0.7	VIOLET	45		G6F7.62Q45D	
F8-90Q-200-D 8mm	 Length: 8mm Radius: 5.1mm Curvature: 90° Chord: 7.21mm Diameter: 200µ	6-0	0.7	VIOLET	30		G6F8Q30D	
		5-0	1	VIOLET	30		G5F8Q30D	
F8.51-90Q-200-D 8.51mm	 Length: 8.51mm Radius: 5.42mm Curvature: 90° Chord: 7.67mm Diameter: 200µ	6-0	0.7	VIOLET	20		G6F8.51Q20D	
		6-0	0.7	VIOLET	45		G6F8.51Q45D	
		6-0	0.7	UNDYED	45		G6F8.51Q45DU	
		5-0	1	UNDYED	45		G5F8.51Q45DU	


7/16 Circle		Micropoint Spatula Needles						Micro Needles
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
F5.5-158V-200-D 5.5mm	 Length: 5.5mm Radius: 2.00mm Curvature: 158° Chord: 3.93mm Diameter: 200µ	7-0	0.5	VIOLET	45		G7F5.5V45D	





3/8 Circle		X-Ultima Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
U6.2-135C-200	 Length: 6.2mm Radius: 2.63mm Curvature: 135° Chord: 4.86mm Diameter: 200µ	10-0	0.2	VIOLET	15	MONOFILAMENT	GAU6.2C15	
6.2mm								




3/8 Circle		X-Trocar Needles					Micro Needles	
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
E12-135C-300	 Length: 12mm Radius: 5.01mm Curvature: 135° Chord: 9.42mm Diameter: 300µ	6-0	0.7	VIOLET	70		G6E12C70	
12mm		5-0	1	VIOLET	70		G5E12C70	


1/2 Circle		Taper Point Plus Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q17U		17mm	6-0	0.7	VIOLET	70		G6Q17U70
			5-0	1	VIOLET	70		G5Q17U70
			5-0	1	UNDYED	70		G5Q17U70U
			5-0	1	VIOLET	70	X-BLACK	G5Q17U70B
			4-0	1.5	VIOLET	70		G4Q17U70
			4-0	1.5	UNDYED	70		G4Q17U70U
			4-0	1.5	VIOLET	4 X 45	CONTROL RELEASE	G4Q17U45
			4-0	1.5	VIOLET	8 X 45	CONTROL RELEASE	G4Q17U45
			4-0	1.5	VIOLET	8 X 70	CONTROL RELEASE	G4Q17U70
			4-0	1.5	VIOLET	70	X-BLACK	G4Q17U70B
			3-0	2	VIOLET	70		G3Q17U70
			3-0	2	UNDYED	70		G3Q17U70U
			3-0	2	VIOLET	8 X 45	CONTROL RELEASE	G3Q17U45
			2-0	3	VIOLET	70		G2Q17U70
Q17U-D		17mm	3-0	2	VIOLET	70	G3Q17U70D	


1/2 Circle		Taper Point Plus Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q22U		22mm	4-0	1.5	VIOLET	70		G4Q22U70
			4-0	1.5	UNDYED	70		G4Q22U70U
			4-0	1.5	VIOLET	70	X-BLACK	G4Q22U70B
			3-0	2	VIOLET	70		G3Q22U70
			3-0	2	UNDYED	70		G3Q22U70U
			3-0	2	VIOLET	4 X 45	CONTROL RELEASE	G3Q22U45
			3-0	2	VIOLET	70	X-BLACK	G3Q22U70B
			2-0	3	VIOLET	70		G2Q22U70
			2-0	3	UNDYED	70		G2Q22U70U
			2-0	3	VIOLET	8 X 45	CONTROL RELEASE	G2Q22U45
			2-0	3	VIOLET	70	X-BLACK	G2Q22U70B




1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R17U		17mm	6-0	0.7	VIOLET	75		G6R17U75
			5-0	1	VIOLET	75		G5R17U75
			5-0	1	UNDYED	75		G5R17U75U
			4-0	1.5	VIOLET	75		G4R17U75
			4-0	1.5	UNDYED	75		G4R17U75U
			4-0	1.5	VIOLET	90		G4R17U90
			4-0	1.5	VIOLET	70	HEAVY NEEDLE	G4R17U70H
			3-0	2	VIOLET	75		G3R17U75
			3-0	2	VIOLET	70	X-BLACK	G3R17U70B
			3-0	2	VIOLET	70	HEAVY NEEDLE	G3R17U70H
			2-0	3	VIOLET	75		G2R17U75
			2-0	3	UNDYED	75		G2R17U75U
			2-0	3	VIOLET	70	HEAVY NEEDLE	G2R17U70H
			0	3.5	VIOLET	70		G0R17U70
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0R17U70H
			1	4	VIOLET	70	HEAVY NEEDLE	G1R17U70H



1/2 Circle		Taper Point Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R17U-D		5-0	1	VIOLET	90		G5R17U90D
		3-0	2	VIOLET	70		G3R17U70D
R18U		4-0	1.5	UNDYED	70		G4R18U70U
		4-0	1.5	VIOLET	75		G4R18U75
R20U		4-0	1.5	VIOLET	75		G4R20U75
		3-0	2	VIOLET	75		G3R20U75
		2-0	3	VIOLET	75		G2R20U75
R22U		5-0	1	VIOLET	70		G5R22U70
		4-0	1.5	VIOLET	75		G4R22U75
		4-0	1.5	UNDYED	75		G4R22U75U
		3-0	2	VIOLET	20		G3R22U20
		3-0	2	VIOLET	75		G3R22U75
		3-0	2	VIOLET	90		G3R22U90
		3-0	2	VIOLET	120		G3R22U120
		3-0	2	VIOLET	70	HEAVY NEEDLE	G3R22U70H
		2-0	3	VIOLET	20		G2R22U20





1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R22U		22mm	2-0	3	VIOLET	45		G2R22U45
			2-0	3	VIOLET	75		G2R22U75
			2-0	3	VIOLET	120		G2R22U120
			2-0	3	VIOLET	70	HEAVY NEEDLE	G2R22U70H
			0	3.5	VIOLET	75		G0R22U75
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0R22U70H
			1	4	VIOLET	70	HEAVY NEEDLE	G1R22U70H
R22U-D		22mm	3-0	2	VIOLET	90		G3R22U90D
			2-0	3	VIOLET	90		G2R22U90D
R26U		26mm	5-0	1	VIOLET	75		G5R26U75
			4-0	1.5	VIOLET	45		G4R26U45
			4-0	1.5	VIOLET	70		G4R26U70
			4-0	1.5	VIOLET	90		G4R26U90
			4-0	1.5	VIOLET	8 X 70	CONTROL RELEASE	G4R26U70
			4-0	1.5	VIOLET	70	X-BLACK	G4R26U70B
			3-0	2	VIOLET	20		G3R26U20


1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R26U		26mm	3-0	2	VIOLET	70		G3R26U70
			3-0	2	UNDYED	70		G3R26U70U
			3-0	2	VIOLET	90		G3R26U90
			3-0	2	UNDYED	90		G3R26U90U
			3-0	2	VIOLET	120		G3R26U120
			3-0	2	VIOLET	8 X 45	CONTROL RELEASE	G3R26U45
			3-0	2	VIOLET	70	X-BLACK	G3R26U70B
			3-0	3	UNDYED	70	X-BLACK	G3R26U70BU
			3-0	3	VIOLET	8 X 45	CONTROL RELEASE	G3R26U45B
			3-0	3	VIOLET	8 X 70	CONTROL RELEASE	G3R26U70B
			3-0	2	VIOLET	70	HEAVY NEEDLE	G3R26U70H
			3-0	2	UNDYED	70	HEAVY NEEDLE	G3R26U70HU
			3-0	2	VIOLET	90	HEAVY NEEDLE	G3R26U90H
			2-0	3	VIOLET	20		G2R26U20
			2-0	3	UNDYED	45		G2R26U45U
			2-0	3	VIOLET	70		G2R26U70


1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R26U		26mm	2-0	3	UNDYED	70		G2R26U70U
			2-0	3	VIOLET	90		G2R26U90
			2-0	3	UNDYED	90		G2R26U90U
			2-0	3	VIOLET	120		G2R26U120
			2-0	3	VIOLET	8 X 45	CONTROL RELEASE	G2R26U45
			2-0	3	VIOLET	70	X-BLACK	G2R26U70B
			2-0	3	VIOLET	45	HEAVY NEEDLE	G2R26U45H
			2-0	3	VIOLET	70	HEAVY NEEDLE	G2R26U70H
			2-0	3	UNDYED	70	HEAVY NEEDLE	G2R26U70HU
			2-0	3	VIOLET	90	HEAVY NEEDLE	G2R26U90H
			0	3.5	VIOLET	70		G0R26U70
			0	3.5	UNDYED	70		G0R26U70U
			0	3.5	VIOLET	90		G0R26U90
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0R26U70H
			0	3.5	UNDYED	70	HEAVY NEEDLE	G0R26U70HU
			0	3.5	VIOLET	90	HEAVY NEEDLE	G0R26U90H


1/2 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R26U 	1	4	VIOLET	90		G1R26U90	
	1	4	VIOLET	70	HEAVY NEEDLE	G1R26U70H	
	1	4	VIOLET	90	HEAVY NEEDLE	G1R26U90H	
R26U-D 	3-0	2	VIOLET	90		G3R26U90D	
	2-0	3	VIOLET	90		G2R26U90D	
	2-0	3	VIOLET	150		G2R26U150D	
	0	3.5	VIOLET	70		G0R26U70D	
	1	4	VIOLET	90		G1R26U90D	
R27U 	3-0	2	VIOLET	75		G3R27U75	
	3-0	2	UNDYED	75		G3R27U75U	
	2-0	3	VIOLET	75		G2R27U75	
	2-0	3	UNDYED	75		G2R27U75U	
	0	3.5	VIOLET	75		G0R27U75	
	0	3.5	VIOLET	90		G0R27U90	
	1	4	VIOLET	70		G1R27U75	
	1	4	UNDYED	75		G1R27U75U	



1/2 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R27U 	1	4	VIOLET	90		G1R27U90	
	1	4	VIOLET	75	HEAVY NEEDLE	G1R27U75H	
	2	5	VIOLET	75		GXR27U75	
R30U 	4-0	1.5	VIOLET	70		G4R30U70	
	3-0	2	VIOLET	75		G3R30U75	
	3-0	2	UNDYED	75		G3R30U75U	
	3-0	2	VIOLET	90		G3R30U90	
	2-0	3	VIOLET	75		G2R30U75	
	2-0	3	UNDYED	75		G2R30U75U	
	2-0	3	VIOLET	90		G2R30U90	
	2-0	3	VIOLET	70	HEAVY NEEDLE	G2R30U70H	
	0	3.5	VIOLET	70		G0R30U70	
	0	3.5	UNDYED	70		G0R30U70U	
	0	3.5	VIOLET	90		G0R30U90	
	0	3.5	VIOLET	70	HEAVY NEEDLE	G0R30U70H	
	1	4	UNDYED	75		G1R30U75U	


1/2 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R30U 	1	4	UNDYED	90		G1R30U90U	
	1	4	VIOLET	70	HEAVY NEEDLE	G1R30U70H	
	2	5	VIOLET	70		GXR30U70	
R30U-D 	2-0	3	VIOLET	70		G2R30U70D	
R31U 	3-0	2	VIOLET	70		G3R31U70	
	2-0	3	VIOLET	70		G2R31U70	
	0	3.5	VIOLET	70		G0R31U70	
R36U 	4-0	1.5	VIOLET	90		G4R36U90	
	3-0	2	VIOLET	70		G3R36U70	
	3-0	2	VIOLET	90		G3R36U90	
	2-0	3	VIOLET	70		G2R36U70	
	2-0	3	VIOLET	90		G2R36U90	
	0	3.5	VIOLET	70		G0R36U70	
	1	4	VIOLET	70		G1R36U70	




1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R37U		37mm	4-0	1.5	VIOLET	75		G4R37U75
			4-0	1.5	VIOLET	90	HEAVY NEEDLE	G4R37U90H
			3-0	2	VIOLET	70		G3R37U70
			3-0	2	UNDYED	70		G3R37U70U
			3-0	2	VIOLET	90		G3R37U90
			3-0	2	UNDYED	90		G3R37U90U
			3-0	2	VIOLET	70	HEAVY NEEDLE	G3R37U70H
			3-0	2	UNDYED	70	HEAVY NEEDLE	G3R37U70HU
			3-0	2	VIOLET	90	HEAVY NEEDLE	G3R37U90H
			3-0	2	UNDYED	90	HEAVY NEEDLE	G3R37U90HU
			2-0	3	VIOLET	45		G2R37U45
			2-0	3	VIOLET	75		G2R37U75
			2-0	3	UNDYED	75		G2R37U75U
			2-0	3	VIOLET	90		G2R37U90
			2-0	3	UNDYED	90		G2R37U90U
			2-0	3	VIOLET	70	HEAVY NEEDLE	G2R37U70H




1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R37U		37mm	2-0	3	UNDYED	70	HEAVY NEEDLE	G2R37U70HU
			2-0	3	VIOLET	75	HEAVY NEEDLE	G2R37U75H
			2-0	3	VIOLET	90	HEAVY NEEDLE	G2R37U90H
			2-0	3	UNDYED	90	HEAVY NEEDLE	G2R37U90HU
			0	3.5	VIOLET	75		G0R37U75
			0	3.5	UNDYED	75		G0R37U75U
			0	3.5	VIOLET	90		G0R37U90
			0	3.5	VIOLET	150	LOOP	G0R37U150L
			0	3.5	VIOLET	30	HEAVY NEEDLE	G0R37U30H
			0	3.5	VIOLET	45	HEAVY NEEDLE	G0R37U45H
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0R37U70H
			0	3.5	UNDYED	70	HEAVY NEEDLE	G0R37U70HU
			0	3.5	VIOLET	75	HEAVY NEEDLE	G0R37U75H
			0	3.5	VIOLET	90	HEAVY NEEDLE	G0R37U90H
			0	3.5	UNDYED	90	HEAVY NEEDLE	G0R37U90HU
1	4	VIOLET	75		G1R37U75			



1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R37U		37mm	1	4	UNDYED	75		G1R37U75U
			1	4	VIOLET	90		G1R37U90
			1	4	UNDYED	90		G1R37U90U
			1	4	VIOLET	70	HEAVY NEEDLE	G1R37U70H
			1	4	UNDYED	70	HEAVY NEEDLE	G1R37U70HU
			1	4	VIOLET	75	HEAVY NEEDLE	G1R37U75H
			1	4	VIOLET	90	HEAVY NEEDLE	G1R37U90H
			1	4	UNDYED	90	HEAVY NEEDLE	G1R37U90HU
			2	5	VIOLET	70		GXR37U70
			2	5	VIOLET	75		GXR37U75
			2	5	VIOLET	90		GXR37U90
			2	5	VIOLET	70	HEAVY NEEDLE	GXR37U70H
			2	5	UNDYED	70	HEAVY NEEDLE	GXR37U70HU
			2	5	VIOLET	75	HEAVY NEEDLE	GXR37U75H
			2	5	VIOLET	90	HEAVY NEEDLE	GXR37U90H
			2	5	UNDYED	90	HEAVY NEEDLE	GXR37U90HU



1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R37U-D	 37mm	1	4	VIOLET	90		G1R37U90D
R40U	 40mm	3-0	2	VIOLET	75		G3R40U75
		3-0	2	VIOLET	90		G3R40U90
		3-0	2	UNDYED	90		G3R40U90U
		3-0	2	VIOLET	70	HEAVY NEEDLE	G3R40U70H
		3-0	2	VIOLET	90	HEAVY NEEDLE	G3R40U90H
		3-0	2	UNDYED	90	HEAVY NEEDLE	G3R40U90HU
		2-0	3	VIOLET	70		G2R40U70
		2-0	3	VIOLET	75		G2R40U75
		2-0	3	UNDYED	75		G2R40U75U
		2-0	3	VIOLET	90		G2R40U90
		2-0	3	UNDYED	90		G2R40U90U
		2-0	3	VIOLET	70	HEAVY NEEDLE	G2R40U70H
		2-0	3	UNDYED	70	HEAVY NEEDLE	G2R40U70HU
		2-0	3	VIOLET	90	HEAVY NEEDLE	G2R40U90H


1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R40U		40mm	2-0	3	UNDYED	90	HEAVY NEEDLE	G2R40U90HU
			0	3.5	VIOLET	75		G0R40U75
			0	3.5	VIOLET	90		G0R40U90
			0	3.5	VIOLET	100		G0R40U100
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0R40U70H
			0	3.5	UNDYED	70	HEAVY NEEDLE	G0R40U70HU
			0	3.5	VIOLET	90	HEAVY NEEDLE	G0R40U90H
			0	3.5	UNDYED	90	HEAVY NEEDLE	G0R40U90HU
			0	3.5	VIOLET	120	HEAVY NEEDLE	G0R40U120H
			0	3.5	VIOLET	150	HEAVY NEEDLE/LOOP	G0R40U150HL
			1	4	VIOLET	75		G1R40U75
			1	4	UNDYED	75		G1R40U75U
			1	4	VIOLET	90		G1R40U90
			1	4	UNDYED	90		G1R40U90U
			1	4	VIOLET	100		G1R40U100
			1	4	VIOLET	70	HEAVY NEEDLE	G1R40U70H



1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R40U		40mm	1	4	VIOLET	90	HEAVY NEEDLE	G1R40U90H
			1	4	UNDYED	90	HEAVY NEEDLE	G1R40U90HU
			1	4	VIOLET	120	HEAVY NEEDLE	G1R40U120H
			1	4	VIOLET	150	HEAVY NEEDLE/LOOP	G1R40U150HL
			2	5	VIOLET	75		GXR40U75
			2	5	VIOLET	90		GXR40U90
			2	5	VIOLET	100		GXR40U100
			2	5	VIOLET	70	HEAVY NEEDLE	GXR40U70H
			2	5	VIOLET	90	HEAVY NEEDLE	GXR40U90H
			2	5	UNDYED	90	HEAVY NEEDLE	GXR40U90HU
R40U-D		40mm	3-0	2	VIOLET	150	LOOP	G3R40U150DL
			3-0	2	VIOLET	140	HEAVY NEEDLE	G3R40U140DH
R43U		43mm	0	3.5	VIOLET	70		G0R43U70
			0	3.5	VIOLET	90		G0R43U90
			1	4	VIOLET	70		G1R43U70
			1	4	VIOLET	90		G1R43U90





1/2 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R43U  43mm	1	4	VIOLET	70	HEAVY NEEDLE	G1R43U70H	
	2	5	VIOLET	70		GXR43U70	
	2	5	VIOLET	90		GXR43U90	
	2	5	VIOLET	70	HEAVY NEEDLE	GXR43U70H	
R45U  45mm	3-0	2	VIOLET	75		G3R45U75	
	3-0	2	UNDYED	75		G3R45U75U	
	3-0	2	VIOLET	120		G3R45U120	
	2-0	3	VIOLET	75		G2R45U75	
R48U  48mm	3-0	2	VIOLET	70		G3R48U70	
	3-0	2	UNDYED	90		G3R48U90U	
	2-0	3	VIOLET	70		G2R48U70	
	2-0	3	VIOLET	90		G2R48U90	
	0	3.5	VIOLET	70		G0R48U70	
	0	3.5	VIOLET	90		G0R48U90	
	0	3.5	UNDYED	90		G0R48U90U	
	1	4	UNDYED	90		G1R48U90U	


1/2 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R48U  48mm		0	3.5	VIOLET	100		G0R48U100
		1	4	VIOLET	70		G1R48U70
		1	4	VIOLET	90		G1R48U90
		1	4	VIOLET	100		G1R48U100
		2	5	VIOLET	70		GXR48U70
		2	5	VIOLET	90		GXR48U90
		2	5	UNDYED	90		GXR48U90U
		2	5	VIOLET	100		GXR48U100
		2	5	VIOLET	120		GXR48U120
R60U  60mm		2-0	3	VIOLET	70		G2R60U70
		0	3.5	VIOLET	70		G0R60U70
		1	4	VIOLET	70		G1R60U70
		1	4	VIOLET	150	LOOP	G1R60U150L
		2	5	VIOLET	70		GXR60U70




1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R65U		65mm	2-0	3	VIOLET	70		G2R65U70
			1	4	VIOLET	70		G1R65U70
			1	4	VIOLET	90		G1R65U90
			1	4	VIOLET	120	LOOP	G1R65U120L
			1	4	VIOLET	150	LOOP	G1R65U150L
			1	4	VIOLET	240	LOOP	G1R65U240L
			2	5	VIOLET	90		GXR65U90
			2	5	UNDYED	100		GXR65U100U
			2	5	VIOLET	150		GXR65U150
			2	5	UNDYED	150	LOOP	GXR65U150UL
2	5	VIOLET	240	LOOP	GXR65U240L			
R70U		70mm	1	4	VIOLET	70		G1R70U70






1/2 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R76U		76mm	3-0	2	VIOLET	70		G3R76U70
			2-0	3	VIOLET	70		G2R76U70
			0	3.5	VIOLET	70		G0R76U70
			1	4	VIOLET	75		G1R76U75
			1	4	VIOLET	90		G1R76U90
			1	4	VIOLET	150	LOOP	G1R76U150L
			2	5	VIOLET	70		GXR76U70
			2	5	VIOLET	150	LOOP	GXR76U150L



3/8 Circle		Taper Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R15C		15mm	4-0	1.5	VIOLET	70		G4R15C70
R17C		17mm	4-0	1.5	VIOLET	75		G4R17C75
			4-0	1.5	VIOLET	75	X-BLACK	G4R17C75B
			3-0	2	VIOLET	75		G3R17C75
			3-0	2	VIOLET	75	X-BLACK	G3R17C75B





3/8 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R18C  18mm	5-0	1	VIOLET	70		G5R18C70	
	4-0	1.5	VIOLET	70		G4R18C70	
	3-0	2	VIOLET	70		G3R18C70	
	3-0	2	UNDYED	70		G3R18C70U	
	2-0	3	VIOLET	70		G2R18C70	
R22C  22mm	4-0	1.5	VIOLET	75		G4R22C75	
	3-0	2	VIOLET	75		G3R22C75	
	2-0	3	VIOLET	75		G2R22C75	
R26C  26mm	3-0	2	VIOLET	70		G3R26C70	
	2-0	3	VIOLET	70		G2R26C70	
	0	3.5	VIOLET	70		G0R26C70	
R32C  32mm	3-0	2	VIOLET	75		G3R32C75	
	3-0	2	VIOLET	90		G3R32C90	
	2-0	3	VIOLET	75		G2R32C75	
	2-0	3	VIOLET	90		G2R32C90	
	0	3.5	VIOLET	75		G0R32C75	



3/8 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R32C		32mm	1	4	VIOLET	70	G1R32C70
			2	5	VIOLET	70	GXR32C70




5/8 Circle		Taper Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R170		17mm	3-0	2	VIOLET	90	G3R17090
R260		26mm	3-0	2	VIOLET	70	G3R26070
			2-0	3	VIOLET	70	G2R26070
			2-0	3	UNDYED	70	G2R26070U
			2-0	3	VIOLET	90	G2R26090
			0	3.5	VIOLET	70	G0R26070
			1	4	VIOLET	70	G1R26070
R260-D		26mm	4-0	1.5	VIOLET	70	G4R26070D
			2-0	3	VIOLET	70	G2R26070D
			2-0	3	VIOLET	90	G2R26090D
			0	3.5	VIOLET	70	G0R26070D






5/8 Circle		Taper Point Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R260-D  26mm	0	3.5	VIOLET	90		G0R26090D	
R270  27mm	3-0	2	VIOLET	75		G3R27075	
	2-0	3	VIOLET	75		G2R27075	
	0	3.5	VIOLET	45	HEAVY NEEDLE	G0R27045H	
	0	3.5	VIOLET	75		G0R27075	
	0	3.5	VIOLET	75	HEAVY NEEDLE	G0R27075H	
	1	4	VIOLET	75		G1R27075	
	1	4	VIOLET	75	HEAVY NEEDLE	G1R27075H	
R270-D  27mm	0	3.5	VIOLET	75		G0R27075D	
R310-D  31mm	3-0	2	VIOLET	90		G3R31090D	
R360  36mm	3-0	2	VIOLET	70		G3R36070	
	2-0	3	VIOLET	70		G2R36070	
	0	3.5	VIOLET	70		G0R36070	





5/8 Circle		Taper Point Needles			◎		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R360	 36mm	1	4	VIOLET	70		G1R36070
R360-D	 36mm	0	3.5	VIOLET	90		G0R36090D
R370	 37mm	2-0	3	VIOLET	75		G2R37075
		0	3.5	VIOLET	75		G0R37075
R400	 40mm	2-0	3	VIOLET	70		G2R40070
		0	3.5	VIOLET	70		G0R40070
		0	3.5	VIOLET	90		G0R40090
		1	4	VIOLET	70		G1R40070
		1	4	VIOLET	90		G1R40090
R1100	 110mm	1	4	VIOLET	75		G1R110075



Hook		Taper Point Needles			◎		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R30J	 30mm	2-0	3	VIOLET	70		G2R30J70
		1	4	VIOLET	70		G1R30J70
		1	4	VIOLET	70	HEAVY NEEDLE	G1R30J70H
R32J	 32mm	1	4	VIOLET	75		G1R32J75
R35J	 35mm	0	3.5	VIOLET	70		G0R35J70
R40J-H	 40mm	1	4	VIOLET	70	HEAVY NEEDLE	G1R40J70H



Ski	Taper Point Needles						
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R19S  19mm	4-0	1.5	VIOLET	20		G4R19S20	
	4-0	1.5	VIOLET	120		G4R19S120	
	3-0	2	VIOLET	120		G3R19S120	
	2-0	3	VIOLET	20		G2R19S20	
	2-0	3	VIOLET	120		G2R19S120	




Ski	Taper Point Needles						
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
R23S  23mm	3-0	2	VIOLET	18		G3R23S18	
	3-0	2	VIOLET	120		G3R23S120	
R26S  26mm	3-0	2	VIOLET	20		G3R26S20	
	3-0	2	VIOLET	40	LOOP	G3R26S40L	
	3-0	2	VIOLET	120		G3R26S120	
	2-0	3	VIOLET	20		G2R26S20	
	2-0	3	VIOLET	120		G2R26S120	
	0	3.5	VIOLET	120		G0R26S120	




Straight		Taper Point Needles			⊙		
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R19Y		4-0	1.5	VIOLET	20		G4R19Y20
		4-0	1.5	VIOLET	40	LOOP	G4R19Y40L
		4-0	1.5	VIOLET	70		G4R19Y70
		3-0	2	VIOLET	20		G3R19Y20
		3-0	2	VIOLET	70		G3R19Y70
R19Y-D	 19mm	4-0	1.5	VIOLET	35		G4R19Y35D
R48Y	 48mm	3-0	2	VIOLET	70		G3R48Y70
R51Y	 51mm	2-0	3	VIOLET	90		G2R51Y90
R65Y	 65mm	2-0	3	VIOLET	90		G2R65Y90



1/2 Circle		Taper Point Plus (Heavier Wire Diameter) 					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Q16U-H		3-0	2	VIOLET	70	HEAVY NEEDLE	G3Q16U70H
		2-0	3	VIOLET	70	HEAVY NEEDLE	G2Q16U70H
		0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q16U70H
		1	4	VIOLET	70	HEAVY NEEDLE	G1Q16U70H
Q22U-H		3-0	2	VIOLET	70	HEAVY NEEDLE	G3Q22U70H
		2-0	3	VIOLET	70	HEAVY NEEDLE	G2Q22U70H
		0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q22U70H
		1	4	VIOLET	70	HEAVY NEEDLE	G1Q22U70H
		2	5	VIOLET	70	HEAVY NEEDLE	GXQ22U70H
Q26U-H		3-0	2	VIOLET	70	HEAVY NEEDLE	G3Q26U70H
		2-0	3	VIOLET	70	HEAVY NEEDLE	G2Q26U70H
		2-0	3	UNDYED	70	HEAVY NEEDLE	G2Q26U70HU
		2-0	3	VIOLET	4 X 45	HEAVY NEEDLE	G2Q26U45H
		0	3	VIOLET	4 X 70	HEAVY NEEDLE	G2Q26U70H
		0	3	UNDYED	70	HEAVY NEEDLE	G2Q26U70HU
		0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q26U70H





1/2 Circle		Taper Point Plus (Heavier Wire Diameter)					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q26U-H  26mm	2-0	3	VIOLET	8 X 45	HEAVY NEEDLE	G2Q26U45H	
	1	4	VIOLET	70	HEAVY NEEDLE	G1Q26U70H	
	2	5	VIOLET	70	HEAVY NEEDLE	GXQ26U70H	
Q27U-H  27mm	1	4	VIOLET	90	HEAVY NEEDLE	G1Q27U90H	
Q31U-H  31mm	0	3.5	VIOLET	75	HEAVY NEEDLE	G0Q31U75H	
	1	4	VIOLET	75	HEAVY NEEDLE	G1Q31U75H	
Q36U-H  36mm	3-0	2	VIOLET	70	HEAVY NEEDLE	G3Q36U70H	
	3-0	2	UNDYED	90	HEAVY NEEDLE	G3Q36U90HU	
	2-0	3	VIOLET	70	HEAVY NEEDLE	G2Q36U70H	
	2-0	3	VIOLET	90	HEAVY NEEDLE	G2Q36U90H	
	2-0	3	UNDYED	90	HEAVY NEEDLE	G2Q36U90HU	
	2-0	3	UNDYED	8 X 70	HEAVY NEEDLE	G2Q36U70H	
	0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q36U70H	
	0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q36U70H	
	0	3.5	VIOLET	90	HEAVY NEEDLE	G0Q36U90H	



1/2 Circle		Taper Point Plus (Heavier Wire Diameter)					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
Q36U-H 	36mm	0	3.5	UNDYED	90	HEAVY NEEDLE	G0Q36U90HU
		0	3.5	VIOLET	8 X 45	HEAVY NEEDLE	G0Q36U45H
		1	4	VIOLET	70	HEAVY NEEDLE	G1Q36U70H
		1	4	VIOLET	90	HEAVY NEEDLE	G1Q36U90H
		1	4	VIOLET	90	HEAVY NEEDLE	G1Q36U90H
		1	4	VIOLET	4 X 70	HEAVY NEEDLE	G1Q36U70H
		2	5	VIOLET	90	HEAVY NEEDLE	GXQ36U90H
		2	5	VIOLET	90	HEAVY NEEDLE	GXQ36U90H
		2	5	UNDYED	90	HEAVY NEEDLE	GXQ36U90HU
Q40U-H 	40mm	3-0	2	VIOLET	70	HEAVY NEEDLE	G3Q40U70H
		2-0	3	VIOLET	70	HEAVY NEEDLE	G2Q40U70H
		2-0	3	VIOLET	8 X 45	HEAVY NEEDLE	G2Q40U45H
		0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q40U70H
		0	3.5	VIOLET	90	HEAVY NEEDLE	G0Q40U90H
		1	4	VIOLET	70	HEAVY NEEDLE	G1Q40U70H
		1	4	VIOLET	90	HEAVY NEEDLE	G1Q40U90H




1/2 Circle		Taper Point Plus (Heavier Wire Diameter)					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q40U-H  40mm	1	4	VIOLET	8 X 45	HEAVY NEEDLE	G1Q40U45H	
	2	5	VIOLET	90	HEAVY NEEDLE	GXQ40U90H	
	2	5	UNDYED	90	HEAVY NEEDLE	GXQ40U90HU	
Q45U-H  45mm	1	4	VIOLET	75	HEAVY NEEDLE	G1Q45U75H	
	2	5	VIOLET	70	HEAVY NEEDLE	GXQ45U70H	
Q48U-H  48mm	3-0	2	VIOLET	70	HEAVY NEEDLE	G3Q48U70H	
	2-0	3	VIOLET	70	HEAVY NEEDLE	G2Q48U70H	
	0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q48U70H	
	0	3.5	VIOLET	90	HEAVY NEEDLE	G0Q48U90H	
	0	3.5	UNDYED	90	HEAVY NEEDLE	G0Q48U90HU	
	0	3.5	VIOLET	8 X 45	HEAVY NEEDLE	G0Q48U45H	
	1	4	VIOLET	70	HEAVY NEEDLE	G1Q48U70H	
	1	4	VIOLET	90	HEAVY NEEDLE	G1Q48U90H	
	1	4	UNDYED	90	HEAVY NEEDLE	G1Q48U90HU	
	1	4	VIOLET	8 X 45	HEAVY NEEDLE	G1Q48U45H	
	2	5	VIOLET	90	HEAVY NEEDLE	GXQ48U90H	






1/2 Circle		Taper Point Plus (Heavier Wire Diameter)						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
Q48U-H		48mm	2	5	UNDYED	90	HEAVY NEEDLE	GXQ48U90HU
			2	5	VIOLET	90	HEAVY NEEDLE	GXQ65U90H
Q65U-H		65mm	2	5	VIOLET	4 X 70	HEAVY NEEDLE	GXQ65U70H
			2	5	VIOLET	70	HEAVY NEEDLE	G2Q70U70H
Q70U-H		70mm	2-0	3	VIOLET	70	HEAVY NEEDLE	G2Q70U70H
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0Q70U70H
			1	4	VIOLET	70	HEAVY NEEDLE	G1Q70U70H


1/2 Circle		X-GUARD Blunt Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
G31U		31mm	3-0	2	VIOLET	75		G3G31U75
			2-0	3	VIOLET	75		G2G31U75
G36U		36mm	2-0	3	VIOLET	75		G2G36U75
			2-0	3	VIOLET	75	HEAVY NEEDLE	G2G36U75H
			0	3.5	VIOLET	75		G0G36U75



1/2 Circle		X-GUARD Blunt Point Needles			◎			
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
G36U-H		36mm	0	3.5	VIOLET	90	HEAVY NEEDLE	G0G36U90H
			1	4	VIOLET	90	HEAVY NEEDLE	G1G36U90H
G40U		40mm	0	3.5	VIOLET	75		G0G40U75
			0	3.5	VIOLET	90	HEAVY NEEDLE	G0G40U90H
			1	4	VIOLET	90	HEAVY NEEDLE	G1G40U90H
G45U		45mm	2	5	VIOLET	90		GXG45U90
G48U		48mm	0	3.5	VIOLET	90		G0G48U90
			1	4	VIOLET	90		G1G48U90







1/2 Circle		Blunt Point Needles			○			
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
B30U		30mm	3-0	2	VIOLET	70		G3B30U70
			2-0	3	VIOLET	70		G2B30U70
B35U		35mm	3-0	2	UNDYED	70		G3B35U70U
			2-0	3	VIOLET	70		G2B35U70
			2-0	3	UNDYED	70		G2B35U70U




1/2 Circle		Blunt Point Needles			○			
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
B35U		35mm	2-0	3	VIOLET	90		G2B35U90
			0	3.5	VIOLET	70		G0B35U70
			0	3.5	VIOLET	90		G0B35U90
			1	4	VIOLET	90		G1B35U90
B37U		37mm	3-0	2	UNDYED	90		G3B37U90U
			2-0	3	UNDYED	75		G2B37U75U
			2-0	3	VIOLET	90		G2B37U90
			2-0	3	UNDYED	90		G2B37U90U
			0	3.5	VIOLET	75		G0B37U75
			0	3.5	VIOLET	90		G0B37U90
			0	3.5	UNDYED	90		G0B37U90U
			1	4	VIOLET	75		G1B37U75
			1	4	VIOLET	90		G1B37U90
			1	4	UNDYED	90		G1B37U90U
B40U		40mm	0	3.5	VIOLET	70		G0B40U70
			0	3.5	VIOLET	90		G0B40U90





1/2 Circle		Blunt Point Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
B40U		40mm	1	4	VIOLET	70		G1B40U70
			1	4	VIOLET	90		G1B40U90
B45U		45mm	1	4	VIOLET	90		G1B45U90
			2	5	VIOLET	90		GXB45U90
B48U		48mm	0	3.5	VIOLET	90		G0B48U90
			1	4	VIOLET	90		G1B48U90
			1	4	UNDYED	90		G1B48U90U
B50U		50mm	2-0	3	VIOLET	90		G2B50U90
			0	3.5	VIOLET	90		G0B50U90
			1	4	VIOLET	90		G1B50U90
			1	4	VIOLET	90	HEAVY NEEDLE	G1B50U90H
			2	5	VIOLET	90	HEAVY NEEDLE	GXB50U90H
B64U		64mm	1	4	VIOLET	90		G1B64U90
			2	5	VIOLET	90		GXB64U90



1/2 Circle		Blunt Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
B65U		65mm	1	4	VIOLET	90	G1B65U90
			2	5	VIOLET	90	GXB65U90


3/8 Circle		Blunt Point Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
B64C		64mm	0	3.5	VIOLET	90	G0B64C90
			1	4	VIOLET	70	G1B64C70
B65C		65mm	1	4	VIOLET	100	G1B65C100


1/2 Circle		Taper Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T17U	 17mm	5-0	1	VIOLET	45		G5T17U45
		4-0	1.5	VIOLET	70		G4T17U70
		4-0	1.5	UNDYED	70		G4T17U70U
		3-0	2	VIOLET	70		G3T17U70
		1	4	VIOLET	70		G1T17U70
T22U	 22mm	4-0	1.5	UNDYED	70		G4T22U70U
		2-0	3	VIOLET	70		G2T22U70
T26U	 26mm	3-0	2	VIOLET	70		G3T26U70
		2-0	3	VIOLET	70		G2T26U70
		2-0	3	VIOLET	6 X 45	CONTROL RELEASE	G2T26U45
		0	3.5	VIOLET	70		G0T26U70
T31U	 31mm	2-0	3	VIOLET	75		G2T31U75
		0	3.5	VIOLET	75		G0T31U75
		1	4	VIOLET	75		G1T31U75
T36U	 36mm	3-0	2	VIOLET	90		G3T36U90
		2-0	3	VIOLET	90		G2T36U90





1/2 Circle		Taper Cutting Needles					
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
T36U  36mm	0	3.5	VIOLET	90		G0T36U90	
	1	4	VIOLET	90		G1T36U90	
	2	5	VIOLET	90		GXT36U90	
T37U  37mm	3-0	2	VIOLET	90		G3T37U90	
	2-0	3	VIOLET	75		G2T37U75	
	2-0	3	VIOLET	90		G2T37U90	
	0	3.5	VIOLET	75		G0T37U75	
	0	3.5	VIOLET	90		G0T37U90	
	1	4	VIOLET	75		G1T37U75	
	1	4	VIOLET	90		G1T37U90	
	1	4	UNDYED	90		G1T37U90U	
	2	5	VIOLET	90		GXT37U90	
T40U  40mm	2-0	3	UNDYED	75		G2T40U75U	
	2-0	3	UNDYED	75	SLIM	G2T40U75SU	
	2-0	3	VIOLET	90		G2T40U90	
	0	3.5	VIOLET	75		G0T40U75	



1/2 Circle		Taper Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T40U	 40mm	0	3.5	VIOLET	90		G0T40U90
T45U	 45mm	1	4	VIOLET	75		G1T45U75
T48U	 48mm	2-0	3	VIOLET	90		G2T48U90
		0	3.5	VIOLET	70		G0T48U70
		0	3.5	VIOLET	90		G0T48U90
		1	4	VIOLET	75		G1T48U75
		1	4	VIOLET	90		G1T48U90
		1	4	VIOLET	150	LOOP	G1T48U150L
		2	5	VIOLET	75		GXT48U75
		2	5	VIOLET	90		GXT48U90
		2	5	VIOLET	150	LOOP	GXT48U150L
		6	8	VIOLET	75		GVIT48U75
		6	8	VIOLET	150		GVIT48U150


3/8 Circle		Taper Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T17C  17mm		5-0	1	UNDYED	75		G5T17C75U
		4-0	1.5	UNDYED	70		G4T17C70U
		3-0	2	UNDYED	75		G3T17C75U
T22C  22mm		2-0	3	VIOLET	20		G2T22C20


5/8 Circle		Taper Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T170  17mm		4-0	1.5	UNDYED	75		G4T17075U




Hook		Taper Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T30J-H  30mm		2-0	3	VIOLET	120	HEAVY NEEDLE	G2T30J120H


Hook		Taper Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
T30J-H	 30mm	0	3.5	VIOLET	70	HEAVY NEEDLE	G0T30J70H
		0	3.5	VIOLET	120	HEAVY NEEDLE	G0T30J120H
		1	4	VIOLET	70	HEAVY NEEDLE	G1T30J70H
T30J-D-H	 30mm	2-0	3	VIOLET	60	HEAVY NEEDLE	G2T30J60DH
		0	3.5	VIOLET	45	HEAVY NEEDLE	G0T30J45DH
		0	3.5	VIOLET	60	HEAVY NEEDLE	G0T30J60DH
T31J	 31mm	1	4	VIOLET	75		G1T31J75




Progressive		Taper Cutting Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
T17P	 17mm	4-0	1.5	VIOLET	45	X-BLACK	G4T17P45	
T24P	 24mm	2-0	3	VIOLET	120		G2T24P120	



1/2 Circle		Reverse Cutting Precision Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC18U	 18mm	5-0	1	UNDYED	45		G5PC18U45U	
		4-0	1.5	UNDYED	45		G4PC18U45U	




3/8 Circle		Reverse Cutting Precision Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC16C	 16mm	6-0	0.7	UNDYED	45		G6PC16C45U	
		5-0	1	UNDYED	45		G5PC16C45U	
		3-0	2	UNDYED	45		G3PC16C45U	




3/8 Circle		Reverse Cutting Precision Point Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
PC18C		18mm	4-0	1.5	UNDYED	70		G4PC18C70U
			3-0	2	UNDYED	70		G3PC18C70U
PC19C		19mm	5-0	1	UNDYED	45		G5PC19C45U
			4-0	1.5	UNDYED	45		G4PC19C45U
			3-0	2	UNDYED	45		G3PC19C45U
PC24C		24mm	3-0	2	UNDYED	45		G3PC24C45U


1/2 Circle		X-Prime Reverse Needles						
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
X16U		16mm	5-0	1	UNDYED	45		G5X16U45U
			4-0	1.5	UNDYED	45		G4X16U45U
			4-0	1.5	UNDYED	45		G4X16U45U


3/8 Circle		X-Prime Reverse Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X16C	 16mm	5-0	1	UNDYED	45		G5X16C45U
		4-0	1.5	UNDYED	45		G4X16C45U
X19C	 19mm	5-0	1	UNDYED	45		G5X19C45U
		4-0	1.5	UNDYED	45		G4X19C45U
		4-0	1.5	UNDYED	45	SLIM	G4X19C45SU
		3-0	2	UNDYED	45		G3X19C45U
X24C	 24mm	4-0	1.5	UNDYED	70		G4X24C70U
		3-0	2	UNDYED	45		G3X24C45U
		2-0	3	UNDYED	45		G2X24C45U
X26C	 26mm	3-0	2	UNDYED	70	PRIMIPASS	G3X26C70U



3/8 Circle		X-Prime Reverse Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
X30C	 30mm	2-0	3	VIOLET	70		G2X30C70




3/8 Circle		Conventional Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K16C	 16mm	6-0	0.7	UNDYED	45		G6K16C45U
		5-0	1	UNDYED	45		G5K16C45U
		4-0	1.5	UNDYED	45		G4K16C45U
		4-0	1.5	VIOLET	75		G4K16C75
K19C	 19mm	4-0	1.5	UNDYED	45		G4K19C45U
		3-0	2	UNDYED	45		G3K19C45U




Straight		Conventional Cutting Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K25Y	 25mm	3-0	2	VIOLET	70		G3K25Y70
K48Y	 48mm	3-0	2	UNDYED	70		G3K48Y70U




Straight		Conventional Cutting Needles			▲		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
K60Y  60mm		4-0	1.5	UNDYED	70		G4K60Y70U
		3-0	2	VIOLET	70		G3K60Y70
		3-0	2	UNDYED	70		G3K60Y70U
		2-0	3	VIOLET	70		G2K60Y70
		2-0	3	UNDYED	70		G2K60Y70U




3/8 Circle		Conventional Cutting Precision Point Needles			▲		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
PK16C  16mm		4-0	1.5	UNDYED	45		G4PK16C45U




3/8 Circle		X-Prime Cutting Needles			▲		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D16C  16mm		5-0	1	UNDYED	45		G5D16C45U
		4-0	1.5	UNDYED	45		G4D16C45U
D19C  19mm		4-0	1.5	UNDYED	75		G4D19C75U
		3-0	2	UNDYED	70		G3D19C70U




3/8 Circle		X-Prime Cutting Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
D26C	 26mm	4-0	1.5	UNDYED	70		G4D26C70U
		4-0	1.5	UNDYED	75		G4D26C75U
		3-0	2	UNDYED	75		G3D26C75U
D36C	 36mm	3-0	2	UNDYED	75		G3D36C75U





1/2 Circle		X-Trocar Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E17U	 17mm	5-0	1	VIOLET	70		G5E17U70
		4-0	1.5	VIOLET	70		G4E17U70
		3-0	2	VIOLET	70		G3E17U70
		2-0	3	VIOLET	70		G2E17U70
		0	3.5	VIOLET	70		G0E17U70
		1	4	VIOLET	70		G1E17U70
E20U-H	 20mm	2-0	3	VIOLET	70	HEAVY NEEDLE	G2E20U70H






1/2 Circle		X-Trocar Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
E22U		22mm	4-0	1.5	VIOLET	70		G4E22U70
			3-0	2	VIOLET	70		G3E22U70
			2-0	3	VIOLET	70		G2E22U70
			0	3.5	VIOLET	70		G0E22U70
			1	4	VIOLET	70		G1E22U70
			1	4	UNDYED	70		G1E22U70U
			E26U		26mm	3-0	2	VIOLET
2-0	3	VIOLET				70		G2E26U70
2-0	3	UNDYED				70		G2E26U70U
2-0	3	VIOLET				70	HEAVY NEEDLE	G2E26U70H
0	3.5	VIOLET				70		G0E26U70
0	3.5	VIOLET				70	HEAVY NEEDLE	G0E26U70H
1	4	VIOLET				70		G1E26U70
1	4	VIOLET				70	HEAVY NEEDLE	G1E26U70H
1	4	UNDYED				70	HEAVY NEEDLE	G1E26U70HU



1/2 Circle		X-Trocar Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
E30U		30mm	3-0	2	VIOLET	70		G3E30U70
			2-0	3	VIOLET	70		G2E30U70
			0	3.5	VIOLET	70		G0E30U70
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0E30U70H
			1	4	VIOLET	70		G1E30U70
			1	4	VIOLET	90	HEAVY NEEDLE	G1E30U90H
E37U		37mm	3-0	2	VIOLET	70	HEAVY NEEDLE	G3E37U70H
			3-0	2	VIOLET	90		G3E37U90
			3-0	2	UNDYED	90		G3E37U90U
			2-0	3	VIOLET	70		G2E37U70
			2-0	3	VIOLET	70	HEAVY NEEDLE	G2E37U70H
			2-0	3	VIOLET	90		G2E37U90
			2-0	3	UNDYED	90		G2E37U90U
			2-0	3	VIOLET	90	HEAVY NEEDLE	G2E37U90H
			2-0	3	VIOLET	150	LOOP	G2E37U150L
			0	3.5	VIOLET	70		G0E37U70




1/2 Circle		X-Trocar Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E37U  37mm		0	3.5	VIOLET	70	HEAVY NEEDLE	G0E37U70H
		0	3.5	VIOLET	90		G0E37U90
		0	3.5	UNDYED	90		G0E37U90U
		0	3.5	VIOLET	90	HEAVY NEEDLE	G0E37U90H
		1	4	VIOLET	70		G1E37U70
		1	4	VIOLET	70	HEAVY NEEDLE	G1E37U70H
		1	4	VIOLET	90		G1E37U90
		1	4	UNDYED	90		G1E37U90U
		1	4	VIOLET	90	HEAVY NEEDLE	G1E37U90H
		2	5	UNDYED	70	HEAVY NEEDLE	GXE37U70HU
		2	5	VIOLET	90		GXE37U90
		2	5	UNDYED	90		GXE37U90U
		2	5	VIOLET	90	HEAVY NEEDLE	GXE37U90H
E40U  40mm		2-0	3	VIOLET	70		G2E40U70
		2-0	3	VIOLET	70	HEAVY NEEDLE	G2E40U70H
		2-0	3	VIOLET	90		G2E40U90




1/2 Circle		X-Trocar Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
E40U		40mm	0	3.5	VIOLET	70		G0E40U70
			0	3.5	VIOLET	70	HEAVY NEEDLE	G0E40U70H
			0	3.5	UNDYED	70	HEAVY NEEDLE	G0E40U70HU
			0	3.5	VIOLET	90		G0E40U90
			1	4	VIOLET	70		G1E40U70
			1	4	VIOLET	70	HEAVY NEEDLE	G1E40U70H
			1	4	UNDYED	70	HEAVY NEEDLE	G1E40U70HU
E43U		43mm	2-0	3	VIOLET	70		G2E43U70
			2-0	3	VIOLET	70	HEAVY NEEDLE	G2E43U70H
			2-0	3	VIOLET	90		G2E43U90
			0	3.5	VIOLET	70		G0E43U70
			0	3.5	VIOLET	150	LOOP	G0E43U150L
			0	3.5	VIOLET	150	HEAVY NEEDLE, LOOP	G0E43U150HL
			1	4	VIOLET	70	HEAVY NEEDLE	G1E43U70H
			1	4	VIOLET	90		G1E43U90
			1	4	VIOLET	150	LOOP	G1E43U150L

1/2 Circle		X-Trocar Needles					
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E43U-H	 43mm	1	4	VIOLET	150	HEAVY NEEDLE, LOOP	G1E43U150HL
		2	5	VIOLET	70	HEAVY NEEDLE	GXE43U70H
		2	5	VIOLET	90	HEAVY NEEDLE	GXE43U90H
E48U	 48mm	2-0	3	VIOLET	70		G2E48U70
		2-0	3	VIOLET	90		G2E48U90
		0	3.5	VIOLET	70		G0E48U70
		0	3.5	VIOLET	90		G0E48U90
		1	4	VIOLET	70		G1E48U70
		1	4	VIOLET	90		G1E48U90
		1	4	VIOLET	150	LOOP	G1E48U150L
		2	5	VIOLET	70		GXE48U70
		2	5	VIOLET	90		GXE48U90
		2	5	VIOLET	150	LOOP	GXE48U150L
E80U	 80mm	1	4	VIOLET	90		G1E80U90

3/8 Circle		X-Trocar Needles						
	Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code	
E18C		18mm	5-0	1	VIOLET	70		G5E18C70
			4-0	1.5	VIOLET	70		G4E18C70
			4-0	1.5	UNDYED	70		G4E18C70U
			3-0	2	VIOLET	70		G3E18C70
E22C		22mm	4-0	1.5	VIOLET	70		G4E22C70
			3-0	2	VIOLET	70		G3E22C70
			2-0	3	VIOLET	70		G2E22C70
E26C		26mm	3-0	2	VIOLET	70		G3E26C70
			2-0	3	VIOLET	70		G2E26C70
E32C		32mm	3-0	2	VIOLET	70		G3E32C70
			2-0	3	VIOLET	70		G2E32C70
			0	3.5	VIOLET	70		G0E32C70
			1	4	VIOLET	70		G1E32C70
			2	5	VIOLET	70		GXE32C70

5/8 Circle		X-Trocar Needles					
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
E170	 17mm	4-0	1.5	UNDYED	70		G4E17070U

1/2 Circle		Taper Point, Reverse Cutting Needles			 		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R37-C24U-H	 37,24mm	2-0	3	VIOLET	150	HEAVY NEEDLE	G2R37-C24U150H

1/2 Circle		Taper Point, Reverse Cutting Needles			 		
Needle		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
R-C37U-H	 37mm	2-0	3	VIOLET	150	HEAVY NEEDLE	G2R-C37U150H

Precut		Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
		4-0	1.5	VIOLET	12 X 45		G4TY1245
		4-0	1.5	UNDYED	12 X 45		G4TY1245U
		4-0	1.5	VIOLET	3 X 45		G4TY345
		4-0	1.5	VIOLET	6 X 45		G4TY645

Precut						
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
	3-0	2	VIOLET	12 X 45		G3TY1245
	3-0	2	UNDYED	12 X 45		G3TY1245U
	3-0	2	VIOLET	12 X 60		G3TY1260
	3-0	2	VIOLET	3 X 45		G3TY345
	3-0	2	VIOLET	3 X 75		G3TY375
	3-0	2	VIOLET	6 X 45		G3TY645
	3-0	2	VIOLET	6 X 75		G3TY675
	2-0	3	VIOLET	12 X 45		G2TY1245
	2-0	3	UNDYED	12 X 45		G2TY1245U
	2-0	3	VIOLET	12 X 60		G2TY1260
	2-0	3	UNDYED	12 X 75		G2TY1275U
	2-0	3	VIOLET	3 X 45		G2TY345
	2-0	3	VIOLET	3 X 75		G2TY375
	2-0	3	VIOLET	6 X 45		G2TY645
	2-0	3	VIOLET	6 X 75		G2TY675
	0	3.5	VIOLET	12 X 45		G0TY1245



Precut						
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
	0	3.5	UNDYED	12 X 45		G0TY1245U
	0	3.5	VIOLET	12 X 60		G0TY1260
	0	3.5	VIOLET	3 X 45		G0TY345
	0	3.5	VIOLET	3 X 75		G0TY375
	0	3.5	VIOLET	6 X 45		G0TY645
	0	3.5	UNDYED	6 X 45		G0TY645U
	0	3.5	VIOLET	6 X 75		G0TY675
	1	4	VIOLET	12 X 45		G1TY1245
	1	4	UNDYED	12 X 45		G1TY1245U
	1	4	VIOLET	12 X 60		G1TY1260
	1	4	VIOLET	3 X 45		G1TY345
	1	4	VIOLET	3 X 75		G1TY375
	1	4	VIOLET	6 X 45		G1TY645
	1	4	UNDYED	6 X 45		G1TY645U
	1	4	VIOLET	6 X 75		G1TY675
	2	5	VIOLET	12 X 45		GXTY1245

Precut						
Needle	Gauge (USP)	Gauge (Metric)	Suture Color	Length (cm)	Other Information	LINXobere Code
	2	5	UNDYED	12 X 45		GXTY1245U
	2	5	VIOLET	3 X 45		GXTY345
	2	5	VIOLET	3 X 75		GXTY375
	2	5	VIOLET	6 X 75		GXTY675